

schülke -†



Now also in softpacks!

Ready-to-use, disinfectant wipes impregnated with alcohol-free cleansing solution.

mikrozid® sensitive wipes

Our plus:

- alcohol-free
- ready-to-use, disinfectant-impregnated wipes
- very fast-acting and effective
- for alcohol-sensitive surfaces
- now also in „premium quality“ in a softpack!

Areas of application

For the disinfection of medical devices [CE 0297] and other environmental surfaces.

Especially suitable for surfaces sensitive to alcohol (e.g. plexiglass and ultrasound transducers).

Instructions for Use

Wipe surfaces thoroughly (ensuring that they are thoroughly wetted) and allow to take effect. When using on particularly large areas, use several wipes if necessary.

Microbiological effectiveness

mikrozid® sensitive wipes is effective against:

- bacteria • fungi • viruses (e.g. hepatitis B, HIV, rotavirus, vaccinia).

Working concentrations / contact times

Effective against	Contact times
Bacteria, fungi (<i>C. albicans</i>) (in accordance with VAH/DGHM Guidelines for shorter contact times)	1 min.
enveloped viruses* (including HIV, hepatitis B and C)	1 min.
Rotaviruses	1 min.
Polyoma SV40	1 min.
Noro-viruses**	30 min.

* in accordance with RKI Recommendation in Bundesgesundheitsblatt [Federal Health Gazette] 01/2004

** Murines Noro-Virus according to EN 14 476

Product data

Composition:

100 g agent solution contain the following active ingredients:
 0.26 g Quaternary ammonium compounds, benzyl-C12-16-alkyldimethyl, chlorides;
 0.26 g Didecyldimethylammonium chloride, 0.26 g Quaternary ammonium compounds, C12-14-alkyl [(ethylphenyl) methyl] dimethyl, chlorides

Chemical / physical data:

Appearance of solution:	clear, colourless liquid
Density (20 °C):	approx 0.998 g/cm ³
Viscosity:	< 15 s. (run-out time acc. to DIN 53 211)
pH:	approx 6.0 – 8.0

mikrozyd® sensitive wipes

Special instructions

Wearing gloves (e.g. of nitril-rubber) is recommended.
Use disinfectant safely. Always read the label information before use.

Delivery form / unit

Size of trading unit	Unit	Art. Nr.
mikrozyd® sensitive wipes		
200 wipes-dispenser	10 tins	109 184
200 wipes-refill	12 refills	109 190
mikrozyd® sensitive wipes premium		
24 wipes softpack	20 packs	165 710
48 wipes softpack	16 packs	165 711

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

Expert opinions and informations

Please visit our website for an overview of all available literature/reports on mikrozyd® sensitive wipes:
www.schuelke.com

For your individual inquiries:
Customer care

Phone: +49 (0)40-521 00-666
Email: info@schuelke.com

The advantages of the mikrozyd® sensitive wipes premium

- handy softpack format
- premium quality wipe
- extra soft
- high coverage



BAUA: N-30345 (PT2)
N- 30346 (PT4)
CH: CHZN 1336



Schülke & Mayr GmbH is certified according to DIN EN ISO 9001, DIN EN ISO 14001 and DIN EN ISO 13485 (Reg.-No. 004567-MP23) and has a validated environmental system in accordance with the Eco Audit Regulation (Reg.-No. DE-150-00003).

Schülke & Mayr GmbH
22840 Norderstedt, Germany
Phone +49 (0) 40 - 521 00 - 0
Fax +49 (0) 40 - 521 00 - 318
www.schuelke.com
mail@schuelke.com

Schülke & Mayr UK Ltd.
1 Jenkin Road
GB-Sheffield S9 1AT
Phone +44 (0) 1142 - 54 35 - 00
Fax +44 (0) 1142 - 54 35 - 01
mail.uk@schuelke.com

Schülke & Mayr Ges.m.b.H
Seidengasse 9
1070 Wien, Austria
Phone +43 (0) 1 - 523 25 01 - 0
Fax +43 (0) 1 - 523 25 01 - 60
office.austria@schuelke.com