

schülke -†



Microbial filtrated, liquid concentrate based on quaternary ammonium compounds and biguanide for disinfecting surfaces in cleanrooms

perform® advanced concentrate QB

At a glance:

- Fast & effective
- Microbial filtrated (0.2 µm)
- Available in a single-use bottle
- Colour & perfume free
- Efficacy tested to European norms

Applications

Disinfection of surfaces and equipment/devices in clean room areas and other areas with strict hygienic requirements.

Instructions for use

To prepare the active solution, fill a tub or bucket with the appropriate amount of water, pour the concentrate into the water and stir gently.

For example:

1 bottle (40 ml) for 8 L water = 0.5 % solution

Wipe surfaces so they are well moistened and observe contact time.

Microbiological effectiveness

perform® advanced concentrate QB is: bactericidal, levurocidal, virucidal (BVDV, Vaccina).

Working concentrations/contact times:

bactericidal	0,5 % (5 ml/l)	5 Min.
levurocidal	0,5 % (5 ml/l)	5 Min.
virucidal (BVDV/Vaccina)	0,5 % (5 ml/l)	60 Min.
	1,0 % (10 ml/l)	5 Min.

Product data

Composition: 100 g contain: 30 g benzyl-C12-18-alkyldimethyl, ammonium, chlorides, 3 g Polyhexamethylene biguanide (monomer: 1,5-bis(trimethylen)-guanylguanidinium monohydrochloride).

Chemical/ physical data:

Appearance:	colourless liquid
Density at 20 °C:	0,99 g/cm ³
pH value:	ca. 4,0

Labelling according to EC Directives

- C: Corrosive.
 N: Dangerous for the environment.
 R21/22: Harmful in contact with skin and if swallowed.
 R34: Causes burns.
 R43: May cause sensitisation by skin contact.
 R51/53: Toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment.
 S26: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
 S35: This material and its container must be disposed of in a safe way.
 S36/37/39: Wear suitable protective clothing, gloves (e.g. made of butyl rubber) and eye/face protection.
 S45: In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).
 S61: Avoid release to the environment. Refer to special instructions/safety data sheets.



Tested according to EN 1276, 1650, 13697.

perform[®] advanced concentrate QB

Special notes

Use disinfectants safely. Always read the label and product information before use.

The practical single-use bottles are easy to open without having to use a device. They are stable and easy to empty. perform[®] advanced concentrate QB features excellent material compatibility.

Forms of delivery/packaging units

Container sizes	Packaging units	Art. No.
40 ml bottle (further packaged as 5-bottle units)	20 x bottles	112701

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

Expert opinions and information

Please visit our website for an overview of all available literature / reports on perform[®] advanced concentrate QB:
www.schuelke.com

For further information please contact your local representative or local offices.



Schülke & Mayr GmbH is certified according to DIN EN ISO 9001, DIN EN ISO 14001 and DIN EN ISO 13485 (Reg.-No. 004567-MP23) and has a validated environmental system in accordance with the Eco Audit Regulation (Reg.-No. DE-150-00003).

Schülke & Mayr GmbH
22840 Norderstedt, Germany
Phone +49 (0) 40 - 521 00 - 0
Fax +49 (0) 40 - 521 00 - 318
www.schuelke.com
mail@schuelke.com

Schülke & Mayr UK Ltd.
1 Jenkin Road
GB-Sheffield S9 1AT
Phone +44 (0) 1142 - 54 35 - 00
Fax +44 (0) 1142 - 54 35 - 01
mail.uk@schuelke.com

Schülke & Mayr Ges.m.b.H
Seidengasse 9
1070 Wien, Austria
Phone +43 (0) 1 - 523 25 01 - 0
Fax +43 (0) 1 - 523 25 01 - 60
office@schuelke.at