



Sterile, alcohol-based active solution
for disinfecting surfaces in cleanrooms

perform[®] sterile alcohol EP

At a glance:

- Fast, effective and residue free
- Sterile (gamma irradiated) and microbial filtrated (0.2 µm)
- Ready to use in a handy, sterile 500ml spray system
- Colour and perfume free
- Efficacy tested to European norms

Applications

Disinfection of surfaces and equipment/devices in clean rooms. Especially for Grade A and B clean areas, according to:

- EU Guide to Good Manufacturing Practice for medicinal products; Annex 1: Manufacture of sterile medicinal products, 2003.
- Manufacturing Chemists' Association (MCA) "Rules and Guidelines for pharmaceutical manufacturers and distributors", 2003; Annex 1, Manufacture of sterile medicinal products.
- FDA Draft Guidance for Industry: Sterile Drug Products Produced by Aseptic Processing, U.S. Department of Health and Human Services, Food and Drug Administration, August 2003.

Instructions for use

Apply perform[®] sterile alcohol EP undiluted to objects and surfaces and allow it to work, ensuring that all areas are moistened completely. Remove any excess preparation with a disposable cloth. perform[®] sterile alcohol EP is also suitable for disinfecting gloves. Do not use to treat surfaces that are particularly sensitive to alcohol, such as acrylic glass.

Microbiological effectiveness

perform[®] sterile alcohol EP is:
bactericidal, fungicidal, virucidal (BVDV, Vaccinia).

Contact times

bactericidal	1 Min.
fungicidal	1 Min.
virucidal (BVDV, Vaccinia)	1 Min.



Tested according to EN 1276, 1650, 13697.

Product data

Composition: 100 g contain: 25 g Ethanol (94 %), 35 g Propan-1-ol.

Chemical / physical data:

Appearance:	clear, colourless liquid
Density at 20 °C:	0,89 g/cm ³
Flash point (DIN 51755):	27 °C

Labelling according to EC Directives

Xi:	Irritant.
R10:	Flammable.
R41:	Risk of serious damage to eyes.
R67:	Vapours may cause drowsiness and dizziness.
S23:	Do not breathe vapour.
S26:	In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
S35:	This material and its container must be disposed of in a safe way.
S39:	Wear eye/face protection.
S51:	Use only in well ventilated areas.

perform[®] sterile alcohol EP

Special notes

Use disinfectants safely. Always read the label and product information before use.

The sterile 500 ml spray system offers the following advantages:

- Especially handy and space-saving.
- Ergonomically shaped spray pump will not slip out of a gloved hand.
- Still sterile after opening.

The bottles are sterile (with radiation indicator) and double-wrapped individually.

The life of a unit taken from an open container corresponds to the use-by date printed on the container with proper handling and storage.

Wearing of protective gloves (e.g., made of butyl rubber) is recommended.

Forms of delivery/packaging units

Container sizes	Packaging units	Art. No.
500 ml spray bottle	10 x 500 ml	121802

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

Expert opinions and information

Please visit our website for an overview of all available literature / reports on perform[®] sterile alcohol EP:

www.schuelke.com

For further information please contact your local representative or local offices.



Schülke & Mayr GmbH is certified according to DIN EN ISO 9001, DIN EN ISO 14001 and DIN EN ISO 13485 (Reg.-No. 004567-MP23) and has a validated environmental system in accordance with the Eco Audit Regulation (Reg.-No. DE-150-00003).



Schülke & Mayr GmbH
22840 Norderstedt, Germany
Phone +49 (0) 40 - 521 00 - 0
Fax +49 (0) 40 - 521 00 - 318
www.schuelke.com
mail@schuelke.com

Schülke & Mayr UK Ltd.
1 Jenkin Road
GB-Sheffield S9 1AT
Phone +44 (0) 1142 - 54 35 - 00
Fax +44 (0) 1142 - 54 35 - 01
mail.uk@schuelke.com

Schülke & Mayr Ges.m.b.H
Seidengasse 9
1070 Wien, Austria
Phone +43 (0) 1 - 523 25 01 - 0
Fax +43 (0) 1 - 523 25 01 - 60
office@schuelke.at