



Wound irrigation solution for cleansing and moistening of chronic skin wounds.

octenilin® wound irrigation solution

Our Plus

- excellent cleaning performance
- especially well-tolerated by the skin and tissues
- for all wound situations
- for repeated and long-term use
- sterile

Application areas

- rapid, effective cleansing of wounds
- for the removal of wound crusts consisting of necrotic tissue, pathogens, biofilm and fibrinous films
- suitable for difficult to access areas, such as fissures and wound pockets
- for moistening and moisturisation of wounds
- for moisturisation of dressings and wound pads
- for gentle removal of difficult to remove, caked dressings/ wound pads
- for combined use with Negative Pressure Wound Therapy (V.A.C.Ult, KCI)

Product Profile

Moisturising property

Due to its low surface tension, octenilin® Wound Irrigation Solution has excellent moisturising properties and thus cleansing properties on tissue surfaces, even difficult to reach surfaces, such as fissures and wound pockets.

Cleansing power to remove biofilms

Studies on the cleansing power against microorganisms in the biofilm revealed a good performance.

In comparison to wound irrigation solutions such as Ringer or isotonic saline solutions commonly used in hospitals, octenilin® Wound Irrigation Solution shows a significantly better cleansing power and, therefore, a better reduction of pathogens.

Conservation

octenilin® Wound Irrigation Solution is conserved as to considerably reduce the risk of germ conveyance into the adjacent surroundings during application. octenilin® Wound Irrigation Solution thus offers, if applied according to the instructions, an effective germ barrier between the wound and the treated person.

Tissue tolerability and biocompatibility

Non-irritating, non-sensitising, free of pain, without tissue toxicity, no impairment granulation and no impairment of the epithelialisation. The good tolerability of Octenidine has been confirmed by clinical experience over many years and by clinical studies on chronic wounds as well.

Instructions for use

Intensively irrigate and clean the wound with octenilin® Wound Irrigation Solution during each change of dressing (octenilin® Wound Irrigation Solution can be warmed up to body temperature immediately before application). Easily removable coverings can already be removed this way. Dressings and wound coverings additionally soaked with octenilin® Wound Irrigation Solution can be used to solvate fibrin coverings. The application should be repeated until all coverings and necroses can be removed and the wound is optically clean. Rinsing afterwards is generally not necessary. The wound can be further treated with octenilin® wound gel, if necessary.

Best before 8 weeks after the first use.

CE 0297



octenilin® wound irrigation solution

Product data

Ingredients:

Aqua valde purificata, Glycerol, Ethylhexylglycerin, Octenidine HCl

Original condition:

Sterile, sealed. Stability at room temperature, see stability information.

Use before:

8 weeks after first opening. The bottle must be closed after each use.

Chemical-physical data

Color	colourless
Density	ca. 1,0 g/cm ³ / 20 °C
Flash point	not applicable
Form	liquid
pH	ca. 5,4 / 20 °C / (undiluted)

Special advice

Use with caked dressings, i. e. dressings that are difficult to remove

If the wound dressing is difficult to remove, soaking it with octenilin® Wound Irrigation Solution is recommended to facilitate the gentle removal of the dressing without traumatising the wound surface.

Side effects

None observed.

Contraindications

octenilin® Wound Irrigation Solution should not be used in patients with known or suspected allergy to one or more of the ingredients. If in doubt consult a physician. octenilin® wound irrigation solution should not be used on hyaline cartilage, the eyes, ears, nose, urinary bladder and in the abdominal cavity! To prevent possible tissue injury, the product must not be injected or applied to tissues with pressure. Adequate drainage from wound cavities must be provided (e.g. for flexible drain tube).

General safety instructions

- Do not use for infusion or injection!
- Not for oral use!
- Only to be applied to the wound externally.
- Only flawless and undamaged container must be used.
- Use only intact and undamaged packs. Protect the product against exposure to direct sunlight.
- To prevent the introduction of bacteria when using octenilin® wound irrigation solution, ensure that the packaging does not come into contact with the wound.
- Medical device - Keep out of the reach and sight of children!

Information for order

Item	Delivery form	Item no.
octenilin wound irrigation 350 ml FL	10/Carton	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

Related Products

- octenilin® wound gel

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining out high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: www.schuelke.com.

For individual questions:

Customer Care

Phone: +49 40 52100-666

E-Mail: info@schuelke.com



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.



schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 (0) 40 - 52100 - 0
Fax +49 (0) 40 - 52100 - 318
www.schuelke.com
mail@schuelke.com

Schülke & Mayr UK Ltd.
1 Jenkin Road
GB-Sheffield S9 1 AT
Phone +44 (0) 1142 - 5435 - 00
Fax +44 (0) 1142 - 5435 - 01
mail.uk@schuelke.com

Schülke & Mayr Ges.m.b.H.
Seidengasse 9
1070 Wien, Austria
Phone +43 (0) 1 - 5232501-0
Fax +43 (0) 1 - 5232501-60
office.austria@schuelke.com