



gigasept® FF (new) is a high-level disinfectant based on succinic dialdehyde for manual reprocessing.

## gigasept® FF (new)

### Our Plus:

- formaldehyde and GDA-free
- broad microbiological efficacy
- outstanding material compatibility
- 14-days standing time

### Application areas

Disinfection of thermolabile and thermostable instruments, e.g. flexible and rigid endoscopes, water-tight ultrasound transducers, anaesthesia accessories, respirator masks, instruments etc.

### Instructions for use

For professional use only.

gigasept® FF (new) is a concentrate and is diluted with cold water to produce the desired application concentration.

Application example: 1 litre of a 6 % working solution corresponds to 940 ml of water and 60 ml of gigasept® FF (new).

Place endoscopes and instruments for reprocessing in the gigasept® FF (new) solution. Ensure complete covering of the surfaces, also of hollow instruments, and allow to take effect.

Rinse all instruments thoroughly after disinfection, using water of at least drinking water quality, or preferably sterile distilled water or fully deionised water, to completely remove residues of the solution. Please follow the reprocessing recommendations of the instrument manufacturers.

Apart from the manual treatment of endoscopes and instruments, gigasept® FF (new) is also suitable for semiautomated machines that operate by a cyclic process at room temperature.

### Standing time

The working solution can be used for 14 days, for the disinfection of previously cleaned instruments only.

National regulations may limit the timespan during which the solution can be used.

### Microbiological efficacy

Efficacy	Concentration	Contact time
<b>bactericidal, yeasticidal, tuberculocidal</b> <i>acc. to VAH</i> EN 13727, EN 14561 EN 13624, EN 14562 EN 14348, EN 14563	5 % (50 ml/l)	15 min.
	4 % (40 ml/l)	30 min.
	2 % (20 ml/l)	60 min.
<b>mycobactericidal</b> <i>acc. to VAH</i> EN 14348, EN 14563	7 % (70 ml/l)	15 min.
	5 % (50 ml/l)	30 min.
	3 % (30 ml/l)	60 min.
<b>enveloped viruses (incl. HIV, HBV, HCV)</b> <i>acc. to DVV-/RKI-Guideline 12/2014</i>	1 % (10 ml/l)	15 min.
<b>virucidal</b> <i>clean conditions EN14476</i>	6 % (60 ml/l)	60 min.
<b>virucidal</b> <i>acc. to DVV-/RKI-Guideline 12/2014</i>	8 % (80 ml/l)	60 min.
<b>sporicidal (C. difficile)</b> <i>clean conditions</i> prEN 17126	6 % (60 ml/l)	8 h
<b>sporicidal (B. subtilis)</b> <i>clean conditions</i> based on EN 14347	8 % (80 ml/l)	16 h



# gigasept® FF (new)

## Product data

### Composition:

100 g solution contain the following active substances:  
93.9 g Reaction product of DMO-THF, ethanol and water.

### Labelling according to Regulation (EC) no. 648/2004:

< 5 % phosphonates, < 5 % anionic surfactants, < 5 % non-ionic surfactants, perfumes, Methyl isothiazolinones.

### Chemical-physical data

Density	ca. 1.01 g/cm <sup>3</sup> / 20 °C
Color	green
Flash point	38.5 °C / Method: DIN 51755 Part 1
Form	liquid
pH	6.3 - 6.6 / 20 °C
Viscosity, dynamic	No data available

## Related products

- gigazyme®
- gigazyme® X-tra

## Special advice

**Always read the label and product information before use.**

### Material compatibility:

metals, rubber, glass, plastics, flexible endoscopes.

### Interactions:

As a general rule, do not mix gigasept® FF (new) with other products. If gigasept® FF (new) comes into contact with application solutions containing amines, there may be precipitation and irreversible discolorations, for example if cleaning agents containing amines are used during precleaning and are not adequately rinsed off.

### Storage:

Do not store for long periods above 25 °C. Do not use preparation after expiry date.

## Information for order

Item	Delivery form	Item no.
gigasept FF new 2 l-bottles	5 / Carton	on request
gigasept FF new 5 l-canister	1 / Canister	on request

Accessories	Item no.
Can key for 5 + 10 l	on request
measuring cup 500 ml	on request
measuring cup 50 ml	on request
Instrumentenwannen 3, 5, 10, 30 l	on request

## Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining out high quality standards.

## Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: [www.schuelke.com](http://www.schuelke.com).

For individual questions:

Customer Sales Service

Phone: +49 (0)40-521 00-666

Email: [info@schuelke.com](mailto:info@schuelke.com)



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.



A company of the Air Liquide-Group

**schülke Headquarters**  
Schülke & Mayr GmbH  
Robert-Koch-Str. 2  
22851 Norderstedt, Germany  
Phone +49 (0) 40 - 52100 - 0  
Fax +49 (0) 40 - 52100 - 318  
[www.schuelke.com](http://www.schuelke.com)  
[mail@schuelke.com](mailto:mail@schuelke.com)