



Sterile, alcohol based ready to use solution for surface disinfection in cleanrooms.

perform® sterile alcohol IPA (WFI)

Our Plus

- fast, effective and residue free
- sterile (X-Ray irradiated) and microbial filtrated (0.2 µm)
- suitable for cleanrooms grade A/B, ISO class 5/6
- easy to use in a handy, sterile 500 ml spray system
- with WFI (Water for Injection)
- colour and perfume free
- efficacy tested according to European norms (EN)

Application areas

Disinfection of surfaces and equipment / devices in cleanrooms, in particular for class A and B cleanroom areas, according to:

- Annex 1 of EU Guide to Good Manufacturing practice: Manufacture of sterile medicinal products, March 2009.
- FDA Guidance for Industry: Sterile drug products produced by aseptic processing - current Good Manufacturing practice; U.S. Department of Health and Human Services, Food and Drug Administration, September 2004.

Instructions for use

Apply perform® sterile alcohol IPA WFI undiluted to objects and surfaces and allow to take effect. Ensure complete wetting. The product dries off quickly. Maximum application rate 50 ml/sq.m. Do not use to treat surfaces that are particularly sensitive to alcohol, such as acrylic glass. perform® sterile alcohol IPA is also available with purified water (VE-water).

Microbiological efficacy

| Efficacy | Concentration | Contact time |
|-----------------------------------------------|---------------|--------------|
| bactericidal <i>EN1276, EN13697</i> | ready-to-use | 1 min. |

| Efficacy | Concentration | Contact time |
|---------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|--------------|
| yeastcidal <i>EN1650, EN13697</i> | ready-to-use | 1 min. |
| virucidal against enveloped viruses <i>in accordance with DVV (German Registered Association for Combating Viral Diseases)//RKI Guideline</i> | ready-to-use | 30 sec. |



perform[®] sterile alcohol IPA (WFI)

Product data

Composition:
100 g solution contains the following active substances: 63,1 g Propan-2-ol (equates 70% (v/v)).

Chemical-physical data

| | |
|--------------------|------------------------------------|
| Color | colourless |
| Density | ca. 0,88 g/cm ³ / 20 °C |
| Flash point | 22 °C / Method : DIN 51755 Part 1 |
| Form | liquid |
| pH | Not applicable |
| Viscosity, dynamic | > 2,5 mPa*s / 20 °C / Raw material |

Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: www.schuelke.com.
For individual questions:
Application Department
Phone: +49 40 52100-666
E-Mail: info@schuelke.com

Special advice

Use disinfectants safely. Always read the label and product information before use.

The sterile 500 ml spray system offers the following advantages:

- convenient size for handling and storage
- ergonomically shaped spray pump for ease of use
- remains sterile after opening

The bottles are sterile (with irradiation indicator) and doublewrapped individually. perform[®] sterile alcohol IPA comprises alcohol of pharmaceutical quality in accordance with the Ph. Eur.

EU-0023656-0001 1-1

Information for order

| Item | Delivery form | Item no. |
|------------------------------------------------------------|---------------|------------|
| perform [®] sterile alcohol IPA WFI 1l bottle | 6/Carton | on request |
| perform [®] sterile alcohol IPA WFI 500 ml bottle | 10/Carton | on request |

These products are not available in every country. For more information please contact our local subsidiary or distributor.

Related Products

- perform[®] sterile alcohol IPA
- perform[®] sterile dry wipes
- perform[®] sterile mix dry wipes
- perform[®] sterile premium dry wipes

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com