



High performance dual action liquid concentrate for cleaning and disinfection of instruments

## gigasept® instru AF

### Our Plus

- very good cleaning performance
- broad spectrum of microbiological efficacy
- shortened contact time with ultrasonic baths

### Application areas

Suitable for manual disinfecting and cleaning of thermostable and thermolabile instruments of all types (excluding flexible endoscopes). Thanks to its excellent cleaning performance especially in an ultrasonic bath and its comprehensive microbiological efficacy, gigasept® instru AF offers all important product features for instrument reprocessing.

### Instructions for use

gigasept® instru AF is a concentrate and is diluted with cold water to produce the desired application concentration.

#### Dosage:

depending on the degree of effectiveness 1.5 % - 4 %.

#### Application example:

1 litre of a 3 % working solution corresponds to 970 ml of water and 30 ml of gigasept® instru AF.

Immerse instruments to be reprocessed in the gigasept® instru AF solution immediately after use. Ensure complete covering of the surfaces, also of hollow instruments, and allow them to soak in the solution. Rinse all instruments thoroughly after cleaning, using water of at least drinking water quality, or preferably sterile distilled water or fully deionised water, to completely remove residues of the cleaning solution.

Do not exceed indicated contact times and concentrations.

Please follow the reprocessing recommendations of the instrument manufacturers. National regulations may require that cleaning and disinfection are carried out in two separate process steps. Not for final disinfection of semicritical and critical medical devices!

#### Standing time:

A standing time of 7 days applies if the working solution is only used for the disinfection of previously cleaned instruments. If the working solution is used for cleaning, national regulations may limit the timespan during which the solution may be used.

#### Ultrasonic cleaning:

gigasept® instru AF is particularly suitable for use in ultrasonic devices. The chemical stability of the product is maintained in this process. Studies have shown that very good cleaning can be achieved with gigasept® instru AF after only a few seconds. Reference reports can be provided on request. Moreover, the stock solution should not be heated above 40 °C.

### Microbiological efficacy

Efficacy	Concentration	Contact time
<b>bactericidal</b> EN13727, EN14561, in accordance with VAH - dirty conditions	3 % (30 ml/l)	15 min.
<b>tuberculocidal</b> EN14348, EN14563, in accordance with VAH - dirty conditions	3 % (30 ml/l)	15 min.
<b>mycobactericidal</b> EN14348, EN14563, in accordance with VAH - dirty conditions	3 % (30 ml/l)	15 min.
<b>fungicidal</b> EN13624, EN14562, in accordance with VAH - dirty conditions	3 % (30 ml/l)	15 min.
<b>virucidal against enveloped viruses</b> in accordance with DVV (German Registered Association for Combating Viral Diseases)//RKI Guideline - dirty conditions	1 % (10 ml/l)	10 min.
<b>Adenovirus</b> in accordance with DVV (German Registered Association for Combating Viral Diseases)//RKI Guideline - clean conditions	3 % (30 ml/l)	60 min.



# schülke →



High performance dual action liquid concentrate  
for cleaning and disinfection of instruments

## gigasept® instru AF

Efficacy	Concentration	Contact time
<b>Polyoma SV40</b> <i>in accordance with DVV (German Registered Association for Combating Viral Diseases)//RKI Guideline - clean conditions</i>	3 % (30 ml/l)	60 min.

CE 0297



## Product data

### Composition:

100 g solution contains the following active substances:  
15.6 g Cocosporylene diamineguanidine diacetate, 35 g  
Phenoxypropanols, 2.5 g Benzalkoniumchloride..

## Special advice

**Always read the label and product information before use.** For professional use only. Any serious incidents related to the product are to be reported to the manufacturer and competent authority. If stock solutions and concentrates of gigasept<sup>®</sup> instru AF come into contact with other products (especially ones containing aldehyde), precipitations/ discolouration/ losses of efficacy are possible. Therefore, any mixture with other products is to be avoided. We therefore recommend rinsing all application aids and tubs thoroughly before using gigasept<sup>®</sup> instru AF for the first time. Keep container tightly closed. Protect from frost, heat and sunlight. Store in the original container. Do not use gigasept<sup>®</sup> instru AF after the expiry date. For all of the related safety information please refer to the Safety Data Sheet. To view the Safety Data Sheet please visit [www.schuelke.com](http://www.schuelke.com).

## Information for order

Item	Delivery form	Item no.
gigasept instru AF -INT- 2 I FL	5/Carton	107413
gigasept instru AF -INT- 5 I KA	1/Canister	107416

## Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

## Expert opinion and information

Please visit our website for an overview of all available literature/reports on all schulke products online at <http://www.schuelke.com/>

For any queries you might have about the product or using it please call 0114 2543500 or E-mail [mail.uk@schuelke.com](mailto:mail.uk@schuelke.com)

Item no.	GTIN no.
107413	4032651074136
144349	-
107416	4032651074167



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

**schülke UK**  
Schülke & Mayr UK Ltd  
Cygnets House, Meadowhall  
Sheffield, S9 1AT  
Phone +44-114-254 35 00  
Fax +44-114-254 35 01  
[mail.uk@schuelke.com](mailto:mail.uk@schuelke.com)

**schülke Headquarters**  
Schülke & Mayr GmbH  
Robert-Koch-Str. 2  
22851 Norderstedt, Germany  
Phone +49 40 - 52100 - 0  
Fax +49 40 - 52100 - 318  
[www.schuelke.com](http://www.schuelke.com)  
[info@schuelke.com](mailto:info@schuelke.com)