



High level disinfectant based on peracetic acid for medical instruments.

**CE** 0297

## gigasept® PAA

### Our Plus

- Aldehyde-free
- Ready to use
- Effective against bacteria, fungi, viruses, spores
- No activation necessary

### Application area

Intended to be used for manual disinfection or semi-automated disinfection of pre-cleaned flexible endoscopes, ophthalmological equipment, surgical instruments made of stainless steel and heat sensitive medical devices.

### Instructions for use

**Use disinfectants safely. Always read the label and product information before use. For professional use only. Read SDS before use.**

gigasept® PAA is a ready to use product. gigasept® PAA is active immediately and stable over 7 days or 50 cycles of disinfection. The peracetic acid (PAA) concentration in solution can be tested using Schülke PAA test strips. In case of visible contamination, the bath should be renewed. Always ensure the pre-cleaning of the endoscope and medical instruments and a complete wetting of the surfaces. Immerse hollow objects in such a way that air escapes completely from the lumen. After the necessary contact time rinse carefully with sterile or filtered water before re-use. Do not exceed the recommended contact time. A located ventilation over the bath is recommended. Cover the bath of gigasept® PAA as often as possible. Do not dilute gigasept® PAA. Do not mix with other products (cleaner, disinfectants, etc.), particularly with aldehydes or derivatives of halogens (e.g. glutaral, sodium hypochlorite, etc.). Do not use after expiry date.

### Microbiological efficacy

	Test method	Contact time
<b>Bacteria, Yeast and Mycobacteria</b>		
bactericidal	suspension test carrier test EN 13727 EN 14561	5 min.
mycobactericidal	suspension test carrier test EN 14348 EN 14563	5 min.
yeastocidal	suspension test carrier test EN 13624 EN 14562	5 min.
<b>Bacterial spores</b>		
Bacillus cereus	suspension test EN 17126	5 min.
Bacillus subtilis	suspension test EN 17126	5 min.
sporicidal	EN 17126	5 min.
sporicidal (C. difficile)	suspension test EN 17126	5 min.
<b>Fungi and Spores</b>		
fungicidal	EN 13624 EN 14562	5 min.
<b>Virus</b>		
virucidal activity	suspension test carrier test EN 14476 EN 17111	5 min.

Product tested under clean conditions.

Note: Description of the test method you will find on our website.



# gigasept® PAA

## Product data

**Composition:** 100 g solution contains the following active substances: 0,12 g Peracetic acid [CAS n°79-21-0].

## Special advice

### Safety:

**WARNING - irritant product- always refer to the product label and the Safety Data Sheet before use.**

gigasept® PAA is not classified as a dangerous formulation by ingestion (Test OECD 401) or irritation to skin (OECD 404).gigasept® PAA does not generate toxic vapour - the average exposure limit of the ingredients is respected.gigasept® PAA is not sensitising (Test OECD 406). Gloves must be worn when handling disinfectants not intended for the skin.Do not dilute gigasept® PAA. Store and use only in well-ventilated area. Do not mix with other products. Do not use after expiry date. Always close canister tightly after removal of product. Canister with venting cap. Do not keep the venting cap sealed. Store in upright position only. Store at ambient temperature in the original canister. Protect from heat and sunlight. Do not store with alkaline products. Never return used product into the original canister. Only recycle completely emptied packaging. Do not reuse empty container. Dispose of product including packaging as hazardous waste in compliance with local and national regulations.

### Compatibility:

gigasept® PAA can be used in semi-automated endoscope reprocessors that are compatible with peracetic acid. Please check the compatibility of endoscope and reprocessors with paracetic acid.Incompatible materials: iron, toolsteel, aluminium, platinum, nickel -/ chromeplated or non ferrous metal alloys parts (e.g. cooper and brass).Compatible materials: (limited to the maximum contact time; obligatory rinse after contact) polycarbonate, polyurethane, silicone, polysulfone, polyamide, polyester, polystyrene.

### Storage and disposal:

Store and use only in a well ventilated place. Store at ambient temperature in the original canister. Always close canister tightly after removal of product. Canister with venting cap. Do not keep the venting cap sealed. Store in upright position only. Protect from heat and sunlight. Do not store with alkaline products. Never return used product into the original canister. Do not reuse empty container. Dispose of the product including packaging as hazardous waste in compliance with local and national regulations

## Information for order

Item	Delivery form	GTIN	Item no.	
gigasept PAA	5 l KA	4 / carton(s)	4260709380094	195751

## Related products

- **disclean® endo PAA**

## Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining out high quality standards.

## Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: [www.schuelke.com](http://www.schuelke.com).

For individual questions:

Application Department

Phone: +49 40 52100-666

E-Mail: [info@schuelke.com](mailto:info@schuelke.com)



## schülke -+

### schülke Headquarters

Schülke & Mayr GmbH  
Robert-Koch-Straße 2  
22851 Norderstedt | Germany  
Phone +49 40 52100582  
[www.schuelke.com](http://www.schuelke.com)  
[info@schuelke.com](mailto:info@schuelke.com)

### Schulke India Pvt. Ltd.

Office # 603, A Wing, Delphi,  
Orchard Avenue, Hiranandani  
Business Park, Powai, Mumbai  
400 076 Maharashtra | India  
Phone 022-61736600  
[www.schuelke.com](http://www.schuelke.com)  
[india.office@schuelke.com](mailto:india.office@schuelke.com)



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.