



High-Level Instrument Disinfectant. Broad-Spectrum. Antimicrobial

gigasept® OPA

Our Plus

- Effective against bacteria, fungi, viruses and spores
- Effective against GTA - resistant strains of Mycobacteria
- Convenient to use as no activation is required
- Fast & Effective
- Excellent material compatibility
- Non-corrosive solution

Application area

For High level disinfection

Instructions for use

- See directions for use before using the product or any kind of application.
- High Level Disinfection:** gigasept® OPA used or reused as a high-level disinfectant for a maximum of 14 days at 25 °C with an immersion time of at least 5 minutes.

Product data

Ortho phthalaldehyde 0.55 % w/v

Special advice

ALWAYS READ THE LABEL AND USE ONLY AS DIRECTED

- gigasept® OPA should not be used to reprocess any urological instruments, which have been used in the treatment of patients with bladder cancer.
- Use only in ventilated areas, wear suitable protective clothing, gloves and eye / face protection.
- In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- May cause sensitization by inhalation and skin contact.

Information for order

Item	Delivery form	Item no.
gigasept® OPA 5L	1	70001125

Environmental information

schulke products are manufactured with advanced, safe and environmentally friendly production processes while always maintaining high quality standards

Expert opinion and information

Please visit our website for an overview of all available literature/ reports on the product: www.schulke.com



schulke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Straße 2
22851 Norderstedt | Germany
Phone +49 40 52100582
www.schuelke.com
info@schuelke.com

Schulke India Pvt. Ltd.
Office # 603, A Wing, Delphi,
Orchard Avenue, Hiranandani
Business Park, Powai, Mumbai
400 076 Maharashtra | India
Phone 022-61736600
www.schuelke.com
india.office@schuelke.com



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.