



For acute wounds and minor burns, as well as sunburn and itching

CE 0297

octenisept® gel

Our Plus

- moistens and protects the wound
- promotes the natural healing process
- painless on application
- helps to alleviate itching
- contains Octenidine as a preservative, capable of inactivating pathogens in the gel to build a protective barrier against penetrating germs
- colour and perfume free

Application area

octenisept® gel is designed for treating acute wounds such as grazes, cuts, minor burns, and sunburn. It aids optimal wound healing by maintaining a moist environment, which is crucial for skin cell division and migration. The gel forms a protective layer against bacteria, inactivating pathogens and preventing infections. This process supports and accelerates natural wound healing. Additionally, octenisept® gel has a cooling effect that alleviates itching. It is colorless and non-burning upon application.

Product Profile

Tissue tolerability and biocompatibility of octenisept® gel Non-irritating, non-sensitising, pain-free, without tissue toxicity, no impairment of granulation or epithelialisation. The good tolerability of Octenidine has been confirmed by clinical experience over many years and by clinical studies on chronic wounds.

Instructions for use

Clean and disinfect the wound immediately after the injury occurs. Then, cover the wound with a thin layer of octenisept® gel. If required, apply a dressing or plaster. Apply the wound gel 1 to 2 times a day until the wound has healed up completely. For hygienic reasons, the tip of the gel tube must not touch the wound or other surfaces.

Product data

Ingredients:

Aqua purificata, Propylene Glycol, Hydroxyethylcellulose, Octenidine HCl.

Chemical-physical data:

Density: ca. 1,01 g/cm³ / 20 °C

Form: viscous

Flash point: > 61 °C / Method : ISO 2719

Color: nearly colourless

Viscosity, dynamic: ca. 5.050 mPa*s

pH: ca. 6,9 / 20 °C

Special advice

Side effects

None observed to date.

Contraindications

Do not use octenisept® gel continuously for longer than 2 weeks without medical supervision. As allergies can never be excluded, octenisept® gel should not be applied if allergies are known or in the case of a suspected allergy to one of the ingredients. If in doubt consult a physician. For external use only. To prevent possible tissue damage, octenisept® gel must not be applied on cartilage, in the ear or nose and must not be allowed to enter the



octenisept® gel

eye.

General safety instructions

- Not for infusion or injection
- Not for internal use
- With acute skin lesions, for external use only.
- For chronic wounds, use a sterile wound gel (octenilin® wound gel).
- Only use unopened and undamaged packages.

Information for order

Item	Delivery form	GTIN	Item no.
octenisept gel -EN/EL/PT- ml TB 20	20 / carton(s)	4032651216116	121611

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining out high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature/reports on all schulke products online at [http:// www.schuelke.com/](http://www.schuelke.com/)

For any queries you might have about the product or using it please call 0114 2543500 or E-mail mail.uk@schuelke.com



schülke -+

Manufacturer

Schülke & Mayr GmbH
Robert-Koch-Straße 2
22851 Norderstedt | Germany
Phone +49 40 52100582
www.schuelke.com
info@schuelke.com

schülke UK

Schulke & Mayr UK Ltd.
1 Jenkin Road
S9 1AT Sheffield | United Kingdom
Phone +44 114 254 3500
<http://www.schulke.co.uk>
mail.uk@schuelke.com



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.