



Low alcohol rapid disinfection of non-invasive medical devices and surfaces.

mikrozid® universal wipes green line

Our Plus

- 100% plastic-free wipes
- sustainable and climate-neutral
- unique low alcohol formulation - blend of two alcohols
- dual action - cleaning and disinfection with one product
- effective against enveloped viruses
- added surfactants for outstanding cleaning
- compatible with a wide range of materials including: leatherette, tablets and displays
- effective against bacteria within one minute contact time (EN16615 4-field test)
- dermatologically tested

Application areas

Alcohol-based rapid cleaning and disinfection of medical devices and surfaces when short contact times are required and where a material friendly disinfectant is needed, e.g.:

- Patient treatment units
- Examination couches
- Operating tables and nearby work surfaces
- Digital displays
- Surfaces of medical equipment
- Keyboards or touch pads of sensitive equipment, e.g. smartphones or tablets

Instructions for use

Thoroughly wipe the surface with the impregnated wipe and allow to take effect. Make sure to wet surfaces completely and keep them wet for the whole exposure time. Make sure that all visible dirt is removed prior to disinfection. When opening the softpack, make sure that hygiene requirements are maintained. Make sure to close the lid tightly after use. Not suitable for the disinfection of invasive medical devices.

schülke guarantees a shelf life of 1 month after the package has been opened.

Microbiological efficacy

Efficacy	Concentration	Contact time
bactericidal EN13727, EN16615 - clean conditions	ready-to-use	1 min.
bactericidal EN13727 - dirty conditions	ready-to-use	15 sec.
tuberculocidal EN14348 - clean conditions	ready-to-use	5 min.
yeastocidal EN16615 - clean conditions	ready-to-use	1 min.
yeastocidal EN13624 - dirty conditions	ready-to-use	1 min.
virucidal against enveloped viruses EN14476 - dirty conditions	ready-to-use	15 sec.
Norovirus EN14476 - clean conditions	ready-to-use	30 sec.
Norovirus EN14476 - dirty conditions	ready-to-use	60 sec.
Polyoma SV40 in accordance with DVV (German Registered Association for Combating Viral Diseases)//RKI Guideline - dirty conditions	ready-to-use	15 sec.
Rotavirus EN14476 - clean conditions	ready-to-use	15 sec.



mikrozid® universal wipes green line

Product data

Composition:

100 g solution contains the following active substances: 17,4 g

Propan-2-ol, 12,6 g Ethanol (94 % w/w)

Labelling according to Regulation (EC) No. 648/2004: < 5%
anionic surfactants

Special advice

Special instructions

For all of the related safety information please refer to the Safety Data Sheet. To view the Safety Data Sheet please visit www.schuelke.com

Use disinfectants safely. Always read the label and product information before use.

The product is highly compatible with metals and synthetics. Please follow carefully any national fire and explosion protection guidelines when using alcohol-based disinfectants. Only treat small surfaces. The wiping range depends on the room temperature and on the structure of the surface to be disinfected. More information available on request.

schülke guarantees, that with appropriate handling, a shelf life of 1 month after opening the packaging.

The following sizes are available:• mikrozid® universal wipes green line (6 x 114 wipes 18 x 20cm)

Information for order

Item	Delivery form	Item no.
mikrozid univ wipes green line 114 ST SP	6/Carton	70003354

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature/reports on all schulke products online at <http://www.schuelke.com>/For any queries you might have about the product or using it please call 01 14 2543500 or E-mail mail.uk@schuelke.com



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.