



## aspirmatic®

Cleaning, disinfection and total care for all dental system units and amalgam separators. Use daily 1-2 times.

### Our Plus

- suitable for all amalgam separators (eg Dürr, Metasys)
- can be used in all treatment units (eg Sironna, KaVo, Castellini, etc.)
- does not foam

### Application areas

aspirmatic® is used for the daily cleaning, disinfection and care of all dental system units and amalgam separators.

### Instructions for use

Use aspirmatic® in a 2% solution at noon and at the end of the day for a cleaning disinfection. Use 2 l in-use solution of aspirmatic® for every dental system unit and use in combination with the schülke dental application aid. To make up a 2 l in-use solution add 40 ml aspirmatic® to 1960 ml cold water. Dilute the product with cold water to get the required application concentration.

Suck 1 1/2 litres of in-use solution through the suction hoses and pour the remaining 1/2 litre in the rinse basin. Leave to soak for a minimum of 60 minutes, or overnight. We recommend to use or aspirmatic® cleaner 1-2 per week instead of the aspirmatic®.

### Microbiological efficacy

Efficacy	Concentration	Contact time
<b>bactericidal</b> in accordance with VAH - dirty conditions	2 % (20 ml/l)	60 min.
<b>bactericidal</b> EN13727, EN14561 - dirty conditions	2 % (20 ml/l)	60 min.
<b>yeastcidal</b> EN13624, EN14562 - dirty conditions	2 % (20 ml/l)	60 min.

Efficacy	Concentration	Contact time
<b>yeastcidal</b> in accordance with VAH - dirty conditions	2 % (20 ml/l)	60 min.
<b>virucidal against enveloped viruses</b> in accordance with DVV (German Registered Association for Combating Viral Diseases)//RKI Guideline	2 % (20 ml/l)	60 min.

### Certificates

- VAH certificate

CE 0279



## Product data

### Composition:

100 g concentrate contains the following active substances:  
7,5 g dioctyl dimethyl ammonium chloride.

Labelling according to Regulation (EC) No 648/2004: < 5 %  
non-ionic surfactants, < 5 % soap.

### Chemical-physical data

Color	blue
Density	ca. 0,99 g/cm <sup>3</sup> / 20 °C
Flash point	Not applicable
Form	liquid
pH	6,5 - 7,5 / 20 °C
Viscosity, dynamic	No data available

## Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product:

[www.schuelke.com](http://www.schuelke.com).

For individual questions:

Application Department

Phone: +49 40 52100-666

E-Mail: [info@schuelke.com](mailto:info@schuelke.com)

## Special advice

### Always read the label and product information before use.

The UVV (accident prevention regulations) require the persons dealing with disinfectants which are not determined for the skin to wear gloves. Do not use after the expiration date. Do not mix the concentrate or the ready-to-use solution of aspirmatic® and aspirmatic® cleaner.

## Information for order

Item	Delivery form	Item no.
aspirmatic 2 I FL	5/Carton	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

## Application aids

Application aids	Item no.
measuring cup 500 ml	136101
measuring cup 50 ml	136102

## Related Products

- aspirmatic® cleaner
- schülke dental application aid

## Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining out high quality standards.



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

**schülke Headquarters**  
Schülke & Mayr GmbH  
Robert-Koch-Str. 2  
22851 Norderstedt, Germany  
Phone +49 40 - 52100 - 0  
Fax +49 40 - 52100 - 318  
[www.schuelke.com](http://www.schuelke.com)  
[info@schuelke.com](mailto:info@schuelke.com)