



Sterile liquid concentrate based on peracetic acid for disinfecting surfaces in cleanrooms.

perform® sterile concentrate PAA

Our Plus

- fast, effective and sporicidal
- sterile (X-Ray irradiated)
- convenient single use bottle
- colour and perfume free
- residue free
- efficacy tested according to European norms (EN)

Application areas

Disinfection of surfaces and equipment/devices in cleanrooms. Especially for Grade A and B clean areas, according to:

- EU Guidelines to Good Manufacturing Practice for medicinal products, Annex 1: Manufacture of sterile medicinal products; March 2009.
- FDA Guidance for Industry: Sterile drug products produced by aseptic processing - current good manufacturing practice, U.S. Department of Health and Human Services; Food and Drug Administration; September 2004.

Instructions for use

Please ensure compatibility of materials and instruments with peracetic acid.

Prepare the working solution by filling drinking quality water into a tub or bucket. Add the concentrate to the water and stir gently.

For example:

1 flask (40 ml) for 4 litre water (water + product) = 1 % solution.

Wipe the pre-cleaned surfaces until they are well wetted and observe required contact time.

Microbiological efficacy

Efficacy	Concentration	Contact time
bactericidal <i>EN1276, EN13697</i>	0.25 % (2.5 ml/l)	5 min.
yeastocidal <i>EN1650, EN13697</i>	0.2 % (2 ml/l)	15 min.
fungicidal <i>EN1650, EN13697</i>	2 % (20 ml/l)	15 min.
virucidal against enveloped viruses <i>in accordance with DVV</i>	0.5 % (5 ml/l)	10 min.
sporicidal <i>EN13704</i>	1 % (10 ml/l)	5 min.



perform[®] sterile concentrate PAA

Product data

Composition:

100 g concentrate contain the following active substances: 6.3 g Peracetic acid. (CAS n°79-21-0)

Other ingredients: Contains Hydrogen peroxide, acetic acid.

Chemical-physical data

Color	colourless
Density	1,123 g/cm ³ / 20 °C
Flash point	> 100 °C
Form	liquid
pH	< 1 / 20 °C
Viscosity, dynamic	1,47 mPa*s / 20 °C

Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: www.schuelke.com. For individual questions:

Application Department

Phone: +49 40 52100-666

E-Mail: info@schuelke.com

Manufacturer:

Bioxal

ZI Sud Secteur A

Route des Varennes

71100 Chalon sur Saône

France

Special advice

Safety: DANGER - Oxidising, corrosive and irritant product, harmful for the environment - always refer to the product label and the Safety Data Sheet before use.

Compatibility:

Good material compatibility with stainless steel, PTFE, PVDF, PEHD, PP, Viton, glass etc.

Contact with non-ferrous metals (brass, copper) should be avoided.

Storage and disposal:

Recommended storage temperature: 5°C to 30°C.

Store and use only in well-ventilated area. Do not mix with other products. Do not use after expiry date. Keep sterile package tightly closed. Do not use if sterile package is open or damaged. Canister with venting cap. Do not keep the venting cap sealed. Store in upright position only. Store at ambient temperature in the original canister. Protect from heat and sunlight. Do not store with alkaline products. Only recycle completely emptied packaging. Do not reuse empty container. Dispose of the product as hazardous waste including packaging in compliance with local and national regulations.

Information for order

Item	Delivery form	Item no.
perform [®] sterile concentrate PAA 3x40 ml bottle	6/Carton	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

Manufacturer
Bioxal
ZI Sud Secteur A
Route des Varennes
71100 Chalon sur Saône
France

schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com