



intended to be used for moistening and decontamination by physical cleansing of the nasal vestibules as well as for supportive treatment of irritated skin underneath the nasal opening.

octenisan® md nasal gel

Our Plus

- decontamination by physical cleansing of the nasal vestibules
- to support treatment of irritated skin underneath the nasal opening
- for moistening
- excellent compatibility


Application areas

The hydrophilic properties of octenisan® md nasal gel ensure efficient moisturising of the nasal vestibules. The moisturising effects of hydrogel promote a favourable physiological milieu which prevents the settlement of bacteria. Furthermore, restoring the physiological functions of protection and cleaning contribute to effective decontamination of e.g. particle, skin residuals and multiresistant pathogens. Irritated skin underneath the nasal opening is also supplied with moisture through the hydrogel and thus can heal optimally. Epithelial cells must be able to proliferate and migrate for the healing process. This is best possible in a moist wound environment.

Instructions for use

Unless otherwise prescribed, octenisan® md nasal gel should be used up to twice a day as needed. octenisan® md nasal gel is for use in the nasal vestibules and on the skin underneath the nasal opening only. Administration in the nasal vestibules: The gel should be applied e.g. with a cotton swab in a sufficient quantity to the frontal area of the nose (nasal vestibules), so that the surface is completely covered. For this purpose, when used in the nasal vestibules, octenisan® md nasal gel should be inserted into the nostril, the nose should be closed by pressing on the side and then massaging between thumb and finger. This ensures an even distribution of the gel. If needed, the nose should be thoroughly cleaned before with a wet cotton swab. Excessive gel is to be removed with a cotton swab. It should be ensured that a moist film

remains. Direct contact with the opening of the tube must be avoided to prevent contamination with germs. Administration underneath the nasal opening: The gel should be applied e.g. with a cotton swab in a sufficient quantity on the irritated skin underneath the nasal opening, so that the surface is completely covered. It should be ensured that a moist film remains.



CE 0297



octenisan[®] md nasal gel

Product data

Ingredients:

Purified water, Propylene Glycol, Hydroxyethylcellulose, Octenidine HCl

Chemical-physical data

Color	nearly colourless
Density	ca. 1,01 g/cm ³ / 20 °C
Flash point	> 61 °C / Method : ISO 2719
Form	viscous
pH	ca. 6,9 / 20 °C
Viscosity, dynamic	ca. 5.050 mPa*s

Expert opinion and information

For an overview of our products containing octenidine, please refer to our website www.schuelke.com.

For individual questions:

Customer Sales Service

Phone: +49 40 52100-666

E-Mail: info@schuelke.com

Special advice

octenisan[®] md nasal gel must not be used with known or suspected allergy to one of the ingredients (propylene glycol, hydroxyethyl cellulose, octenidine HCl). If in doubt, a physician should be consulted.

The product should not be used in combination with cleansing soaps, ointments, oils, enzymes or similar substances, as this may affect the preservation of the product. It should not be used in combination with PVP iodine, as this can lead to discoloration and impaired antiseptic effect of the PVP iodine.

Do not use the octenisan[®] md nasal gel without medical supervision for longer than 2 weeks without interruption. octenisan[®] md nasal gel must not be used on cartilage, in the eyes or in the ear to avoid possible tissue damage. Do not insert too deep in the nose. Do not use for infusion or injection! Do not swallow. Do not use if tube is damaged. octenisan[®] md nasal gel should be stored between 5 °C and 25 °C and protected from direct sunlight.

Information for order

Item	Delivery form	Item no.
octenisan md Nasengel INT SM 6 ml TB	20/Carton	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining out high quality standards.



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

Manufacturer

Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com

Distributor

Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com

Schülke & Mayr Ges.m.b.H.
Seidengasse 9
1070 Wien, Austria
Phone +43 (0) 1 - 5232501-0
Fax +43 (0) 1 - 5232501-60
www.schuelke.at
office.austria@schuelke.com