



## gigasept® AF

gigasept® AF is an aldehyde-free preparation for cleaning and disinfection of medical devices (e.g. flexible and rigid endoscopes) and breathing masks.

### Our Plus

- contains no aldehydes
- suitable for flexible and rigid endoscopes
- effective against bacteria, yeasts, enveloped viruses
- fresh odor

### Application areas

gigasept® AF is for the manual cleaning and disinfection of surgical and thermolabile instruments such as flexible and rigid endoscopes and breathing masks.

### Instructions for use

For professional use only.

gigasept® AF is a concentrate and is diluted with cold water to produce the desired application concentration. **Dosage:** depending on the degree of effectiveness 1% - 4%.

**Application example:** 1 litre of a 4 % working solution corresponds to 960 ml of water and 40 ml of gigasept® AF. Place endoscopes and instruments for reprocessing in the gigasept® AF solution immediately after use. Ensure complete covering of the surfaces, also of hollow instruments, and allow to take effect. Rinse all instruments thoroughly after cleaning/ disinfection, using water of at least drinking water quality, or preferably sterile distilled water or fully deionised water, to completely remove residues of the solution. Please follow the reprocessing recommendations of the instrument manufacturers. National regulations may require that cleaning and disinfection are carried out in two separate process steps. Do not mix gigasept® AF with other disinfectants and cleaners. Not for final disinfection of semicritical and critical medical devices!

**Standing time:** The working solution can be used for 7 days, when used for the disinfection of previously cleaned instruments only. If the working solution is used for cleaning, national regulations may limit the timespan during which the solution may be used.

### Microbiological efficacy

Efficacy	Concentration	Contact time
<b>bactericidal</b> <i>in accordance with VAH - dirty conditions</i>	1.5 % (15 ml/l)	60 min.
<b>bactericidal</b> <i>in accordance with VAH - dirty conditions</i>	3 % (30 ml/l)	30 min.
<b>bactericidal</b> <i>in accordance with VAH - dirty conditions</i>	4 % (40 ml/l)	15 min.
<b>yeastcidal</b> <i>in accordance with VAH - dirty conditions</i>	1.5 % (15 ml/l)	60 min.
<b>yeastcidal</b> <i>in accordance with VAH - dirty conditions</i>	3 % (30 ml/l)	30 min.
<b>yeastcidal</b> <i>in accordance with VAH - dirty conditions</i>	4 % (40 ml/l)	15 min.
<b>virucidal against enveloped viruses</b> <i>in accordance with DVV (German Registered Association for Combating Viral Diseases)/RKI Guideline - dirty conditions</i>	1 % (10 ml/l)	15 min.

CE 0297



## Product data

### Composition:

100 g solution contains the following active substances the following active ingredients: 15 g Didecyldimethylammoniumchloride, 10 g Phenoxypyrrols, 6.9 g Amino alkyl glycine (2.1 g Amines, n-C10-16-Alkyltrimethylenedi-, reaction products with Chloroacetic acid, 4.8 g Betaines, C12-14-Alkyldimethyl).

### Labelling according to regulation (EC) No. 648/2004:

15 % - 30 % non-ionic surfactants, perfumes

### Chemical-physical data

Color	green
Density	ca. 1,00 g/cm <sup>3</sup> / 20 °C
Flash point	45 °C / Method : DIN 51755 Part 1
Form	liquid
pH	ca. 9,0 / 20 °C
Viscosity, dynamic	not determined

## Special advice

### Ultrasonic bath:

gigasept® AF can be used in ultrasonic bath. The concentration of use under high organic load (*bactericidal, levurocidal, enveloped viruses incl. HIV, HBV, HCV*) amounts 4 % at a contact time of 10 min.

### Security advice:

**Always read the label and product information before use.**

Do not use this preparation after the expiration date lapses.

With the use of 4 % solutions, irritation of the mucosa can occur with accidental contact to the eyes.

### Interactions:

Danger of precipitation / discolouration / loss of efficacy if use solution or concentrate of gigasept® AF in to contact with other products or household cleaners (e.g. aldehyde-based products).

## Information for order

Item	Delivery form	Item no.
gigasept AF 2 l FL	5/ Carton	on request
gigasept AF 5 l KA	1/ Canister	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

## Application aids

Application aids	Item no.
Can key for 5 + 10 l	135810
Decanter funnel	117901
schülke dosing feeder 5 l / 10 l (20ml / stroke)	117101
schülke tap for 5 l / 10 l can	135501

## Related Products

- bath-system 10l
- bath-system 30l
- bath-system 3l
- bath-system 5l
- gigasept® PAA concentrate

## Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining out high quality standards.

## Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: <http://www.schuelke.com/>. For individual questions: Customer Sales Service Phone: +49 40 52100-666 E-Mail: [info@schuelke.com](mailto:info@schuelke.com)



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

**schülke Headquarters**  
Schülke & Mayr GmbH  
Robert-Koch-Str. 2  
22851 Norderstedt, Germany  
Phone +49 40 - 52100 - 0  
Fax +49 40 - 52100 - 318  
[www.schuelke.com](http://www.schuelke.com)  
[info@schuelke.com](mailto:info@schuelke.com)