# schülke +

## mikrocount® combi

Keimindikatoren  
Agar dip slides

the plus of pure  
performance

## Unser Plus:

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## Our plus:

- Evaluation of total germ count
- Simple application
- Fast, safe and easy
- Precise and reliable
- Economical to use
- Separate evaluation of bacteria, yeasts and moulds

## Applications

- Metalworking and related industry
- Paint industry
- Paper and coatings industry
- Textile industry
- Cosmetic industry
- Water treatment
- Applicable for tests according to DIN 55 997\*

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\* Solvents for paints and varnishes - Deionized water - Requirements and methods of test

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## Quality assurance with mikrocount® combi

In addition to production hygiene measures, hygiene guidelines and quality assurance concepts require routine hygiene controls during the production process and documentation of the results. mikrocount® combi provides every operation with individual means of rapid and reliable hygiene controls. These dip slides can be used for testing raw materials, for in-process controls during the production process, for quality control of finished products, and in watermixed coolants.

mikrocount® combi dip slides enable simple sampling and evaluation of the results, even by personnel without any microbiological training. Special laboratory equipment is not necessary.

## The product

mikrocount® combi is a plastic slide coated on one side with TTC-agar (bacterial growth) and on the other side with Rose-bengal-agar (yeast and mould growth).

mikrocount® combi is manufactured under aseptic conditions. The agar surfaces contain no toxic ingredients that could be transferred to the products being tested.

Furthermore, this dip slide is supplied with a transparent, shatter-proof protective tube. Subsequent contamination is therefore prevented, and the results can be evaluated with the tube sealed, thus avoiding the risk of contamination to employees.

With strict production standards, we guarantee that mikrocount® combi is a product of consistent high quality.

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## Product benefits

### Easy to use

Without any preparation time, the slide can be used to test surfaces, semi-solid materials and liquids.

### Quick and convenient

Always handy and reliable to use – even at last minute, e.g. when testing needs to be carried out during night shift.

No special laboratory equipment is necessary.

### Economical

The previous time-consuming laboratory work for preparing nutrient media, producing dilution series and counting colonies is eliminated. In addition, the considerable costs of external laboratory are saved.

### Accurate and reliable

The strict standards to which mikrocount® combi is produced guarantee a product of consistent excellent quality.

### Detection of bacteria and moulds

With mikrocount® combi, bacteria, yeasts and moulds can be separately detected with a single dip slide.

### Simple to use

Sampling and evaluation can be carried out even by personnel who have no training in microbiology.

### Safe

Incubated samples are evaluated in sealed tubes which risks of contamination are avoided.

## Storage and stability

The mikrocount® combi tubes should be stored unopened at room temperature (about + 20 °C) and protected from draught and light. The slides must not be allowed to freeze. The expiration date is marked on the box. Any unused slide showing microbial growth should be discarded. Temperature fluctuations may result in condensation settling at the bottom of the dip slide tube.

This sterile liquid can simply be disposed of down the sink.

The performance and the results of mikrocount® combi are not affected by this, provided that the media are not visually dehydrated.

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## Instructions

1. Unscrew the lid of the container and remove the slide without touching the agar surfaces.
2. For direct testing in circulation systems, dip the slide into the liquid to be tested at an accessible point, or hold the slide in the liquid stream. If the liquid sprays at high pressure, make sure the agar does not detach from the slide. If samples are being taken using a beaker, mix the liquid thoroughly before dipping in the slide into the liquid. The slide should be held in contact with the sample for about 5 – 10 seconds. Both sides of the slide must be completely wetted.  
  
Taking samples by dipping the slide into the liquid shall not be done for high viscous liquids. For high viscous fluids e.g. like paints and also at points that are difficult to access, a sterile swab shall be used to take the sample, which is then streaked out onto the agar surface. The high viscous liquid shall be streaked out evenly as a very thin layer onto the agar.  
  
For surfaces and solid media, press each side of the slide on the surface.
3. Place the slide back in the tube and screw on the top.
4. After incubation for 24 – 48 hours at 27 – 30°C, the TTC agar result (bacteria) is available. In 2 – 4 days the result is usually readable at room temperature. For slow-growing microorganisms, the result should be checked again after 48 hours. Yeasts and moulds grow after an incubation period of 72 hours. In 4 – 7 days the result is usually readable at room temperature.
5. After incubation, compare the colony density on the agar surface with the evaluation chart. The results should be evaluated with the tube sealed. Ideally, the incubation temperature should be as close as possible to practical conditions.

## Disposal

Disposal of used mikrocount® combi is best carried out by burning or autoclaving. We shall be pleased to send you appropriate disposal recommendations.

mikrocount® combi	Art.-Nr.
1 Box à 10 Dip-slides	184 901

mikrocount® combi	Art.-No.
1 box contains 10 dip slides	184 901

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**Packungsbeilage beachten!**

## Interpretation of results

The evaluation chart shows colonies formed on mikrocount® combi which correspond to different degrees of microbial contamination. The figures are shown in cfu\*/ml.

1. Colony growth of incubated samples on both sides of the slide is compared to the pictures in the evaluation chart. The picture that most closely resembles the colony density on the agar slide is taken as the test result.
2. **Bacterial growth**  
Most bacteria grow to form red colonies. Colourless colonies may also form. Caution! These colourless colonies must also be taken into consideration.
3. **Moulds and yeasts**  
It is possible for moulds or yeasts or a mixture of moulds and yeasts to grow on the Rose-bengal-agar.

\* cfu = colony forming units

**Pay attention to the package insert!**

Total Bacterial Count Agar						
Bakterien (KBE/ml)						
Bacteria (cfu/ml)						
per ml	10 <sup>2</sup>	10 <sup>3</sup>	10 <sup>4</sup>	10 <sup>5</sup>	10 <sup>6</sup>	10 <sup>7</sup>

Rose-Bengal-Agar					
Hefen (KBE/ml)					
Yeasts (cfu/ml)					
per ml	10 <sup>2</sup>	10 <sup>3</sup>	10 <sup>4</sup>	10 <sup>5</sup>	10 <sup>6</sup>

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For incubation of mikrocount® combi samples we recommend the small "cultura®" incubator

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## product benefits of cultura®:

- compact enough to use almost anywhere
- easy temperature adjustment
- designed for use with mikrocount® combi or mikrocount® TPC

## cultura® – the versatile small incubator

The cultura® incubator is compact and versatile enough for almost any laboratory or manufacturing setting. The built in tray has room to hold up to 18 mikrocount® combi or mikrocount® TPC dip slide samples. A transparent door allows for viewing of the contents without removing samples from the incubator.

The adjustable temperature is pre-set by the manufacturer to maintain 30 °C, an optimum temperature for incubating mikrocount® dip slides.

Results for bacteria are available after 24 to 48 hours.

The detection of yeast and mold takes slightly longer (72 hours).

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Our recommendations regarding our products are based on in-depth scientific testing in our Research Department; they are given in good faith, but no liability can be derived from them. It is the responsibility of the final product manufacturer to assure that claims made for the final product are in conformance with all applicable local laws. In other respect our Conditions of Sale and Supply apply.

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22840 Norderstedt | Germany  
Phone +49 40 521 00-0 | Fax +49 40 521 00-244  
www.schuelke.com | sai@schuelke.com

A company of the  
Air Liquide Group

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Schülke & Mayr Benelux B.V.
Schülke & Mayr AG
Schülke & Mayr UK Ltd.
### Schülke & Mayr GmbH