

schülke -t



Sterile ready-to-use solution based on non-ionic surfactants for cleaning surfaces in cleanrooms

perform[®] sterile cleaner ND

At a glance:

- Ready to use
- Sterile (gamma irradiated) and microbial filtrated (0.2 µm)
- Available in convenient sterile 500 ml spray system
- Colour and perfume free
- Excellent cleaning performance

Applications

Cleaning of surfaces and equipment/devices in clean rooms. Especially for Grade A and B clean areas, according to:

- EU Guide to Good Manufacturing Practice for medicinal products; Annex 1: Manufacture of sterile medicinal products, 2003.
- Manufacturing Chemists' Association (MCA) "Rules and Guidelines for pharmaceutical manufacturers and distributors", 2003; Annex 1, Manufacture of sterile medicinal products.
- FDA Draft Guidance for Industry: Sterile Drug Products Produced by Aseptic Processing, U.S. Department of Health and Human Services, Food and Drug Administration, August 2003.

Instructions for use

Apply perform[®] sterile cleaner ND undiluted to objects and surfaces and wipe if necessary.

Product data

Labelling according to Regulation (EC) No. 648/2004: < 5 % non-ionic surfactants.

Chemical / physical data:

Appearance:	clear, colourless liquid
Density at 20°C:	1,0 g/cm ³
pH:	ca. 5,0

Special notes

The sterile 500 ml spray system offers the following advantages:

- Convenient and space saving.
- Ergonomically shaped spray pump will not slip out of a gloved hand.
- Still sterile after opening.

The bottles are sterile (with radiation indicator) and double-wrapped individually.

perform[®] sterile cleaner ND

Forms of delivery/packaging units

Container sizes	Packaging units	Art. No.
500 ml spray bottle	10 x 500 ml	182002

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

Expert opinions and information

Please visit our website for an overview of all available literature / reports on perform[®] sterile cleaner ND:
www.schuelke.com

For further information please contact your local representative or local offices.



Schülke & Mayr GmbH is certified according to DIN EN ISO 9001, DIN EN ISO 14001 and DIN EN ISO 13485 (Reg.-No. 004567-MP23) and has a validated environmental system in accordance with the Eco Audit Regulation (Reg.-No. DE-150-00003).



Schülke & Mayr GmbH
22840 Norderstedt, Germany
Phone +49 (0) 40 - 521 00 - 0
Fax +49 (0) 40 - 521 00 - 318
www.schuelke.com
mail@schuelke.com

Schülke & Mayr UK Ltd.
1 Jenkin Road
GB-Sheffield S9 1AT
Phone +44 (0) 1142 - 54 35 - 00
Fax +44 (0) 1142 - 54 35 - 01
mail.uk@schuelke.com

Schülke & Mayr Ges.m.b.H
Seidengasse 9
1070 Wien, Austria
Phone +43 (0) 1 - 523 25 01 - 0
Fax +43 (0) 1 - 523 25 01 - 60
office@schuelke.at