



For moistening, cleansing and decontamination of the nasal vestibules

octenisan® nasal gel

Our Plus

- decontamination by physical cleaning of the nasal vestibules
- for moistening
- excellent compatibility

Application areas

octenisan® Nasal Gel is used for moistening, cleansing and decontamination by physical cleaning of the nasal vestibules. External influences can cause the nose to rapidly dry out and become irritated. This can be remedied by supplying the nasal vestibules with sufficient moisture.



Instructions for use

Using a cotton bud, introduce a sufficient amount into the front part of the nose (nasal vestibules). Introduce octenisan® Nasal Gel into each nostril, close the nose by pressing together the sides of the nasal wings and then massaging between the thumb and forefinger. This ensures that the gel is evenly distributed. If necessary, the nose should be thoroughly cleaned beforehand with a thoroughly wetted cotton bud. Remove excess gel using a swab. Ensure that a film of moisture remains behind. For nasal dryness apply max. twice daily on demand. For decontamination apply octenisan® Nasal Gel twice daily for five days. For use in the nasal vestibules. Take care not to introduce octenisan® Nasal Gel too deep into the nose.



octenisan[®] nasal gel

Product data

Ingredients:

aqua, propylene glycol, hydroxyethylcellulose, octenidine HCl.

Chemical-physical data

Color	nearly colourless
Density	ca. 1,01 g/cm ³ / 20 °C
Flash point	> 61 °C / Method: ISO 2719
Form	viscous
pH	ca. 6,9 / 20 °C
Viscosity, dynamic	ca. 5.050 mPa*s

Special advice

In patients with a known hypersensitivity to one of the ingredients, octenisan[®] Nasal Gel must not be used. Do not use octenisan[®] Nasal Gel in combination with other cleaning gels, ointments, soaps, oils, enzymes or other anionic surfactants. Do not use in combination with PVP-iodine, as this may lead to discolorations and may impair the antiseptic action of PVP-iodine. For external use only.

Only flawless and undamaged container must be used.

Store the tube at room temperature protected from direct sunlight.

Information for order

Item	Delivery form	Item no.
octenisan nasal gel 5ml tube	20/Cartron	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining out high quality standards.

Expert opinion and information

For an overview of our products containing octenidine, please refer to our website www.schuelke.com.

For individual questions:

Customer Care

Phone: +49 40 52100-666

E-Mail: info@schuelke.com



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.



schülke Headquarters

Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 (0) 40 - 52100 - 0
Fax +49 (0) 40 - 52100 - 318
www.schuelke.com
mail@schuelke.com

Schülke & Mayr UK Ltd.
1 Jenkin Road
GB-Sheffield S9 1 AT
Phone +44 (0) 1142 - 5435 - 00
Fax +44 (0) 1142 - 5435 - 01
mail.uk@schuelke.com

Schülke & Mayr Ges.m.b.H.
Seidengasse 9
1070 Wien, Austria
Phone +43 (0) 1 - 5232501-0
Fax +43 (0) 1 - 5232501-60
office.austria@schuelke.com



A company of the Air Liquide-Group