



Microbial filtrated, alcohol-based ready to use solution for disinfecting surfaces in cleanrooms.

perform® advanced alcohol EP

Our Plus

- fast, effective and residue free
- microbial filtrated (0.2 µm)
- ready to use
- suitable for cleanrooms grade C/D
- colour and perfume free
- efficacy tested according to European norms (EN)

Application areas

Disinfection of surfaces and equipment/devices in clean room areas and other areas with strict hygienic requirements.

Instructions for use

Apply perform® advanced alcohol EP undiluted to objects and surfaces and allow to take effect. Ensure complete wetting. The product dries off quickly. perform® advanced alcohol EP is also suitable for disinfecting gloves. Maximum application rate 50 ml/sq.m. Do not use to treat surfaces that are particularly sensitive to alcohol, such as acrylic glass.

Microbiological efficacy

Efficacy	Concentration	Contact time
bactericidal <i>EN1276, EN13697</i>	ready-to-use	1 min.
fungicidal <i>EN1650, EN13697</i>	ready-to-use	2 min.
virucidal against enveloped viruses <i>in accordance with DVV (German Registered Association for Combating Viral Diseases)//RKI Guideline</i>	ready-to-use	30 sec.

Efficacy	Concentration	Contact time
limited spectrum virucidal activity <i>EN14476</i>	ready-to-use	30 sec.



Product data

Composition:

100 g solution contains the following active substances: 25 g Ethanol (94 %), 35 g Propan-1-ol.

Chemical-physical data

Color	colourless
Density	approx. 0,89 g/cm ³ / 20 °C
Flash point	27 °C / Method : DIN 51755 Part 1
Form	liquid
pH	approx. 7 / 20 °C

Special advice

Use disinfectants safely. Always read the label and product information before use.

Information for order

Item	Delivery form	Item no.
perform [®] advanced alcohol EP 1l bottle	10/Carton	on request
perform [®] advanced alcohol EP 10l canister	1/Canister	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

Application aids

Application aids	Item no.
Can key for 5 + 10 l	135810
Decanter funnel	117901
schülke tap for 5 l / 10 l can	135501
Spray head for 1l bottle	180212

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: www.schuelke.com.

For individual questions:

Application Department

Phone: +49 40 52100-666

E-Mail: info@schuelke.com



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com