



Sterile, cleanroom compatible dry cloths.

## perform® sterile dry wipes

### Our Plus

- low-particle wipes
- suitable for cleanrooms A/B
- sterile (X-Ray irradiated)
- flexible use for all products of the perform® range
- double bagged in convenient packs of 10

### Application areas

For disinfection and cleaning of surfaces in areas where low-particle wipes are required. Suitable for cleanrooms A/B, ISO class 5.

### Instructions for use

Apply disinfectant/cleaner on the wipe, then wipe objects and surfaces and allow to take effect. Ensure complete wetting. The compatibility of the perform® sterile dry wipes was tested with all sterile disinfectants of the perform® range. Special attention was given to the chemical stability of the active ingredients of the disinfectants in the presence of the wipe material.

The following products of the perform® range are compatible:

- perform® sterile alcohol EP
- perform® sterile alcohol IPA
- perform® sterile PAA
- perform® sterile concentrate PAA
- perform® sterile concentrate OXY
- perform® sterile concentrate Q2

For detailed information on test methods and parameters (e.g. contact times) please contact us.



# perform<sup>®</sup> sterile dry wipes

## Product data

**Material:** polyester, knitted monofilament  
**Size:** 30 x 30 cm  
**Thickness:** 0,62 - 0,65 mm  
**Edges:** laser cut  
**Surface mass:** 130 g/m<sup>2</sup>  
**Fluid intake:** capacity: 635 %  
**Sterility:** by  $\gamma$ -irradiation

## Information for order

Item	Delivery form	Item no.
perform <sup>®</sup> sterile dry wipes 10 pieces / sachet	10/Carton	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

## Related Products

- perform<sup>®</sup> sterile alcohol EP
- perform<sup>®</sup> sterile alcohol IPA
- perform<sup>®</sup> sterile alcohol IPA (WFI)
- perform<sup>®</sup> sterile cleaner ND
- perform<sup>®</sup> sterile concentrate OXY
- perform<sup>®</sup> sterile concentrate PAA
- perform<sup>®</sup> sterile concentrate Q2
- perform<sup>®</sup> sterile mix dry wipes
- perform<sup>®</sup> sterile PAA
- perform<sup>®</sup> sterile premium dry wipes

## Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

## Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: [www.schuelke.com](http://www.schuelke.com).

For individual questions:

Application Department

Phone: +49 40 52100-666

E-Mail: [info@schuelke.com](mailto:info@schuelke.com)



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

**schülke Headquarters**  
Schülke & Mayr GmbH  
Robert-Koch-Str. 2  
22851 Norderstedt, Germany  
Phone +49 40 - 52100 - 0  
Fax +49 40 - 52100 - 318  
[www.schuelke.com](http://www.schuelke.com)  
[info@schuelke.com](mailto:info@schuelke.com)