



For the application of hand disinfectants, washing and care lotions for 500 ml and 1,000 ml disposable bottles.

## Accessories - Dispenser Type KHK / KHL / KHN touchless

### Drip tray round\*



Catch tray frame hangs over the dispenser housing



500 ml	1000 ml	touchless
60000211	60000215	60000219

### Drip tray\*



Catch tray frame with removeable Plastic Tray, hangs over the dispenser housing



500 ml	1000 ml	touchless
669240	669390	669560

### Drip tray SH 1\*



Catch tray frame with removable plastic tray, for wall mounting, short side towards the wall



500 ml	1000 ml	touchless
	60000204	

### Drip tray SH 3\*



Catch tray frame with removable plastic tray, for wall mounting, long side towards the wall



500 ml	1000 ml	touchless
	669310	

### Locking plate



To protect liquids and dispenser against theft

500 ml	1000 ml	touchless 500 ml / 1000 ml
669620	669720	60000337 / 669555



# Accessories - Dispenser

## Type KHK / KHL / KHN touchless

### Stainless steel pump



Stainless steel pump, different dosage adjustments (ca. 0,75/ 1,0/ 1,2/ 1,5 ml), suitable for liquid soaps, alcoholic hand disinfectants, lotions etc.



500 ml	1000 ml	touchless 500 ml / 1000 ml
669601	669701	60000230 60000229

### Plastic pump



Plastic pump with straight suction tube, not autoclavable, different dosage adjustments (ca. 0,7/ 1,0/ 1,2 ml), suitable for liquid soaps, alcoholic hand disinfectants, lotions etc.



500 ml	1000 ml	touchless
669611	669711	-

### Single-use pump



Plastic disposable hygiene pump, hygienic packaged in sealed pouch



500 ml	1000 ml	touchless
60000233	60000232	-

### Angle plate



Angle plate for fixing dispenser parallel to the wall, Aluminium

500 ml	1000 ml und touchless
669319	669320

### Bracket system for Dräger rails



Stainless steel adapter plate for dispenser with aluminium housing

500 ml	1000 ml und touchless
669313	669314

### Bracket system for round tubular frames



Stainless steel adapter plate for dispenser with aluminium housing, bracket for round tubular frame 20 – 35 mm diameter

500 ml	1000 ml und touchless
669311	669312

### Signal frame



Signal frame for dispenser, hung over fitted dispenser backplate

Alternative colours upon request!

### Explanation



autoclavable if supplied with stainless steel pump

134° C, 2 bar



dishwater safe

60° C



single use



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 par 1 German Drug Law and Certificates of GMP-Zertifikaten for medicinal products

**schülke Headquater**  
Schülke & Mayr GmbH  
Robert-Koch-Str. 2  
22851 Norderstedt  
Germany  
Phone +49 40 52100-0  
Fax +49 40 52100-318  
www.schuelke.com  
mail@schuelke.com