

**kodan® tincture forte colourless**

**No Change Service!**

Version  
01.00

Revision Date:  
06.09.2016

Date of last issue: -  
Date of first issue: 06.09.2016

**SECTION 1: Identification of the substance/mixture and of the company/undertaking**

**1.1 Product identifier**

Trade name : kodan® tincture forte colourless

**1.2 Relevant identified uses of the substance or mixture and uses advised against**

Use of the Sub-  
stance/Mixture : Disinfectants

**1.3 Details of the supplier of the safety data sheet**

Manufacturer/ Supplier : Schülke & Mayr GmbH  
Robert-Koch-Str. 2  
  
22851 Norderstedt  
Germany  
Telephone: +49 (0)40/ 52100-0  
Telefax: +49 (0)40/ 52100318  
mail@schuelke.com  
www.schuelke.com

E-mail address of person  
responsible for the  
SDS/Contact person : Application Department HI  
+49 (0)40/ 521 00 8800  
ADHI@schuelke.com  
(Schülke & Mayr UK Ltd.: +44-1142543500)

**1.4 Emergency telephone number**

Emergency telephone num-  
ber : UK Poisons Emergency number: 0870 600 6266

**SECTION 2: Hazards identification**

**2.1 Classification of the substance or mixture**

**Classification (REGULATION (EC) No 1272/2008)**

|  |  |
|--|--|
| Flammable liquids, Category 3                                | H226: Flammable liquid and vapour.       |
| Eye irritation, Category 2                                   | H319: Causes serious eye irritation.     |
| Specific target organ toxicity - single exposure, Category 3 | H336: May cause drowsiness or dizziness. |

**2.2 Label elements**

Hazard pictograms :



**2.3 Other hazards**

This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

Vapours are heavier than air and may spread along floors.

**SECTION 3: Composition/information on ingredients**

**kodan® tincture forte colourless****No Change Service!**Version  
01.00Revision Date:  
06.09.2016Date of last issue: -  
Date of first issue: 06.09.2016**3.2 Mixtures**

Chemical nature : Solution of the following substances with harmless additives.

**Hazardous components**

| Chemical name | Index-Number<br>CAS-No.<br>EC-No.<br>Registration number      | Classification  | Concentration<br>(% w/w) |
|---------------|---|---|--------------------------|
| Propan-2-ol   | 603-117-00-0<br>67-63-0<br>200-661-7<br>01-2119457558-25-XXXX | Flam. Liq. 2; H225<br>Eye Irrit. 2; H319<br>STOT SE 3; H336   | 45                       |
| Propan-1-ol   | 603-003-00-0<br>71-23-8<br>200-746-9<br>01-2119486761-29-XXXX | Flam. Liq. 2; H225<br>Eye Dam. 1; H318<br>STOT SE 3; H336   | 10                       |
| Biphenyl-2-ol | 604-020-00-6<br>90-43-7<br>201-993-5                          | Eye Irrit. 2; H319<br>Eye Dam. 1; H318<br>STOT SE 3; H335<br>Aquatic Acute 1;<br>H400<br>Aquatic Chronic 1;<br>H410 | < 0,5                    |

For explanation of abbreviations see section 16.

**SECTION 4: First aid measures****4.1 Description of first aid measures**

General advice : Take off all contaminated clothing immediately.  
 If inhaled : Move to fresh air. If symptoms persist, call a physician.  
 In case of eye contact : In the case of contact with eyes, rinse immediately with plenty of water and seek medical advice.  
 If swallowed : Do NOT induce vomiting. Drink water as a precaution. Consult a physician if necessary.

**4.2 Most important symptoms and effects, both acute and delayed**

Symptoms : Treat symptomatically.,

**4.3 Indication of any immediate medical attention and special treatment needed**

Treatment : For specialist advice physicians should contact the Poisons Information Service.

**SECTION 5: Firefighting measures****5.1 Extinguishing media**

Suitable extinguishing media : Dry powder, Carbon dioxide (CO<sub>2</sub>), Water spray jet, Alcohol-resistant foam  
 Unsuitable extinguishing media : High volume water jet

**kodan® tincture forte colourless**

**No Change Service!**

Version                      Revision Date:  
01.00                        06.09.2016

Date of last issue: -  
Date of first issue: 06.09.2016

---

**5.2 Special hazards arising from the substance or mixture**

Specific hazards during fire-fighting                      : Cool closed containers exposed to fire with water spray.

Specific risk from the substance or the product itself, its combustion products or evolved gases                      : Vapours may form explosive mixtures with air.

**5.3 Advice for firefighters**

Special protective equipment for firefighters                      : In the event of fire, wear self-contained breathing apparatus.

**SECTION 6: Accidental release measures**

**6.1 Personal precautions, protective equipment and emergency procedures**

Personal precautions                      : Ensure adequate ventilation. Remove all sources of ignition.

**6.2 Environmental precautions**

Environmental precautions                      : Avoid subsoil penetration.

**6.3 Methods and material for containment and cleaning up**

Methods for cleaning up                      : Wipe up with absorbent material (e.g. cloth, fleece).  
Soak up with inert absorbent material (e.g. sand, silica gel, acid binder, universal binder, sawdust).

**6.4 Reference to other sections**

See chapter 13

**SECTION 7: Handling and storage**

**7.1 Precautions for safe handling**

Advice on safe handling                      : Do not spray on a naked flame or any incandescent material.  
Keep away from sources of ignition - No smoking. Keep away from children.

Advice on protection against fire and explosion                      : The hot product gives off combustible vapours. Keep away from sources of ignition - No smoking.

Hygiene measures                      : Keep away from food and drink.

**7.2 Conditions for safe storage, including any incompatibilities**

Requirements for storage areas and containers                      : Store at room temperature in the original container. Keep at temperature not exceeding 25 °C.

Further information on storage conditions                      : Keep away from direct sunlight. Keep container tightly closed.

Advice on common storage                      : Do not store together with oxidising agents.

**7.3 Specific end use(s)**

Specific use(s)                      : none

**SECTION 8: Exposure controls/personal protection**

---

**kodan® tincture forte colourless****No Change Service!**Version  
01.00Revision Date:  
06.09.2016Date of last issue: -  
Date of first issue: 06.09.2016**8.1 Control parameters****Occupational Exposure Limits**

| Components  | CAS-No. | Value type (Form of exposure) | Control parameters                 | Basis    |
|-------------|---------|-------------------------------|------------------------------------|----------|
| Propan-2-ol | 67-63-0 | Permissible exposure limit    | 200 ppm<br>500 mg/m <sup>3</sup>   | TRGS 900 |
|             |         | Ceiling Limit Value           | 400 ppm<br>1.000 mg/m <sup>3</sup> | TRGS 900 |
| Propan-1-ol | 71-23-8 | Permissible exposure limit    | 400 ppm<br>980 mg/m <sup>3</sup>   | OSHA     |
|             |         | Permissible exposure limit    | 200 ppm<br>500 mg/m <sup>3</sup>   | OSHA     |

**Derived No Effect Level (DNEL) according to Regulation (EC) No. 1907/2006:**

| Substance name | End Use | Exposure routes | Potential health effects              | Value                  |
|----------------|---------|-----------------|---------------------------------------|------------------------|
| Propan-2-ol    | Workers | Skin contact    | Long-term exposure, Systemic effects  | 888 mg/kg              |
|                | Workers | Inhalation      | Long-term exposure, Systemic effects  | 500 mg/m <sup>3</sup>  |
| Propan-1-ol    | Workers | Skin contact    | Long-term exposure, Systemic effects  | 136 mg/kg              |
|                | Workers | Inhalation      | Long-term exposure, Systemic effects  | 268 mg/m <sup>3</sup>  |
|                | Workers | Inhalation      | Short-term exposure, Systemic effects | 1723 mg/m <sup>3</sup> |

**Predicted No Effect Concentration (PNEC) according to Regulation (EC) No. 1907/2006:**

| Substance name | Environmental Compartment               | Value      |
|----------------|---|------------|
| Propan-2-ol    | Fresh water                             | 140,9 mg/l |
|                | Marine water                            | 140,9 mg/l |
|                | Fresh water sediment                    | 552 mg/kg  |
|                | Marine sediment                         | 552 mg/kg  |
|                | Soil                                    | 28 mg/kg   |
| Propan-1-ol    | Fresh water                             | 10 mg/l    |
|                | Marine water                            | 1 mg/l     |
|                | Soil                                    | 2,2 mg/kg  |
|                | Marine sediment                         | 2,28 mg/kg |
|                | Fresh water sediment                    | 22,8 mg/kg |
|                | Effects on waste water treatment plants | 96 mg/l    |
|                | Intermittent use/release                | 10 mg/l    |

**8.2 Exposure controls**

**kodan® tincture forte colourless**

**No Change Service!**

Version                      Revision Date:  
01.00                        06.09.2016

Date of last issue: -  
Date of first issue: 06.09.2016

---

**Personal protective equipment**

- Eye protection                      : If splashes are likely to occur, wear:  
Safety glasses with side-shields conforming to EN166
- Protective measures                : Avoid contact with eyes.

**SECTION 9: Physical and chemical properties**

**9.1 Information on basic physical and chemical properties**

- Appearance                        : liquid  
Colour                                : brown  
Odour                                 : alcohol-like  
Odour Threshold                  : not determined  
pH                                      : Not applicable  
Melting point/freezing point    : < -5 °C  
Decomposition temperature      : No data available  
Boiling point/boiling range     : ca. 80 °C  
Flash point                         : 24 °C, DIN 51755 Part 1  
Evaporation rate                  : No data available  
Flammability (solid, gas)        : Not applicable  
Upper explosion limit            : Propan-2-ol: 12 %(V)  
Lower explosion limit            : Propan-2-ol: 2 %(V)  
Vapour pressure                  : ca. 42 hPa, 20 °C  
Relative vapour density         : No data available  
Density                              : ca. 0,90 g/cm<sup>3</sup>, 20 °C  
Solubility(ies)  
    Water solubility                : in all proportions , 20 °C  
Partition coefficient: n-  
octanol/water                      : Not applicable  
Auto-ignition temperature      : Propan-2-ol: 425 °C  
Flow time                          : < 15 s at 20 °C, DIN 53211  
Explosive properties              : No data available  
Oxidizing properties              : No data available

**9.2 Other information**

No data available

**SECTION 10: Stability and reactivity**

**10.1 Reactivity**

No dangerous reaction known under conditions of normal use.

**10.2 Chemical stability**

The product is chemically stable.

**10.3 Possibility of hazardous reactions**

Vapours may form explosive mixture with air.

**10.4 Conditions to avoid**

Heat, flames and sparks.

**10.5 Incompatible materials**

Strong oxidizing agents,

**10.6 Hazardous decomposition products**

**kodan® tincture forte colourless**

**No Change Service!**

Version                      Revision Date:  
01.00                        06.09.2016

Date of last issue: -  
Date of first issue: 06.09.2016

---

None reasonably foreseeable.

**SECTION 11: Toxicological information**

**11.1 Information on toxicological effects**

**Acute toxicity**

**Product:**

Acute oral toxicity                      : Acute toxicity estimate: > 15.000 mg/kg  
Acute inhalation toxicity                : Acute toxicity estimate: 46,3 mg/l  
Acute dermal toxicity                     : Acute toxicity estimate: > 5.000 mg/kg

**Skin corrosion/irritation**

**Product:**

Patch Test 24 Hrs. on humans: No detectable irritation., Influence on wound healing (rat): No negative influence was observed at any time during the healing process

**Serious eye damage/eye irritation**

**Product:**

Rabbit, Causes serious eye irritation.

**Respiratory or skin sensitisation**

**Components:**

**Propan-2-ol:**

Did not cause sensitisation on laboratory animals. Buehler Test, Guinea pig

**Propan-1-ol:**

Does not cause skin sensitisation. Guinea pig, Maximisation Test

**Biphenyl-2-ol:**

Did not cause sensitisation on laboratory animals. Maximisation Test, Guinea pig

**Germ cell mutagenicity**

**Components:**

**Propan-2-ol:**

Germ cell mutagenicity- Assessment                      : Animal testing did not show any mutagenic effects.

**Propan-1-ol:**

Germ cell mutagenicity- Assessment                      : Not mutagenic in Ames Test

**Biphenyl-2-ol:**

Germ cell mutagenicity- Assessment                      : Not mutagenic in Ames Test

**Carcinogenicity**

**Components:**

**Propan-2-ol:**

Carcinogenicity - Assessment                                : Animal testing did not show any carcinogenic effects.

**Propan-1-ol:**

Carcinogenicity - Assessment                                : Animal testing did not show any carcinogenic effects.

**Biphenyl-2-ol:**

**kodan® tincture forte colourless**

**No Change Service!**

Version 01.00  
Revision Date: 06.09.2016

Date of last issue: -  
Date of first issue: 06.09.2016

Carcinogenicity - Assessment : No data available

**Reproductive toxicity**

**Components:**

**Propan-2-ol:**

Reproductive toxicity - Assessment : Animal testing did not show any effects on fertility.

**Propan-1-ol:**

Effects on fertility : Rat, Inhalation, NOAEL: 8,6 mg/l

Effects on foetal development : Rat, Inhalation, NOAEL: 8,6 mg/l

Reproductive toxicity - Assessment : Animal testing did not show any effects on fertility.

**Biphenyl-2-ol:**

Reproductive toxicity - Assessment : No data available

**STOT - single exposure**

**Product:**

May cause drowsiness or dizziness.

**STOT - repeated exposure**

**Components:**

**Propan-2-ol:**

Based on available data, the classification criteria are not met.

**Propan-1-ol:**

The substance or mixture is not classified as specific target organ toxicant, repeated exposure.

**Biphenyl-2-ol:**

No data available

**Aspiration toxicity**

No data available

**SECTION 12: Ecological information**

**12.1 Toxicity**

**Components:**

**Propan-2-ol:**

Toxicity to fish : LC50 (Leuciscus idus): > 100 mg/l, 48 h, static test, Raw material, literature value

Toxicity to daphnia and other aquatic invertebrates : EC50 (Daphnia magna): > 100 mg/l, 48 h, static test, Raw material, literature value

Toxicity to algae : EC50 (Desmodesmus subspicatus (green algae)): > 100 mg/l, 72 h, static test, Raw material, literature value

**Propan-1-ol:**

Toxicity to fish : LC50 (Fish): 3.200 mg/l, 96 h

Toxicity to daphnia and other aquatic invertebrates : EC50 (Daphnia magna (Water flea)): 3.642 mg/l, 48 h

Toxicity to algae : NOEC (Chlorella pyrenoidosa (aglae)): 1.150 mg/l, 48 h

Toxicity to daphnia and other aquatic invertebrates (Chronic toxicity) : NOEC: > 100 mg/l, 21 d, Daphnia magna (Water flea), OECD Test Guideline 211

**kodan® tincture forte colourless**

**No Change Service!**

Version 01.00  
Revision Date: 06.09.2016

Date of last issue: -  
Date of first issue: 06.09.2016

**Biphenyl-2-ol:**

Toxicity to fish : LC50 (Pimephales promelas (fathead minnow)): 5,99 mg/l, 96 h  
 Toxicity to daphnia and other aquatic invertebrates : EC50 (Daphnia magna): 1,5 mg/l, 24 h  
 Toxicity to algae : EC50 (Desmodesmus subspicatus (green algae)): 0,98 mg/l, 72 h  
 M-Factor (Acute aquatic toxicity) : 1  
 Toxicity to daphnia and other aquatic invertebrates (Chronic toxicity) : NOEC: 0,009 mg/l , 21 d, Daphnia magna (Water flea)  
 M-Factor (Chronic aquatic toxicity) : 1

**12.2 Persistence and degradability**

**Product:**

Biodegradability : Readily biodegradable., OECD 301D / EEC 84/449 C6

**Components:**

**Propan-2-ol:**

Biodegradability : Readily biodegradable.

**Propan-1-ol:**

Biodegradability : Readily biodegradable.

**Biphenyl-2-ol:**

Biodegradability : Biodegradation: > 75 %, Exposure time: 28 d, OECD 301B/ ISO 9439/ EEC 84/449 C5

**12.3 Bioaccumulative potential**

**Components:**

**Propan-2-ol:**

Bioaccumulation : No bioaccumulation is to be expected (log Pow <= 4).

Partition coefficient: n-octanol/water : log Pow: 0,05 (20 °C), OECD Test Guideline 107

**Propan-1-ol:**

Bioaccumulation : Bioaccumulation is unlikely.

Partition coefficient: n-octanol/water : log Pow: 0,43

**Biphenyl-2-ol:**

Bioaccumulation : Bioconcentration factor (BCF): 21,07, Bioaccumulation is unlikely.

Partition coefficient: n-octanol/water : log Pow: 3,18

**12.4 Mobility in soil**

**Components:**

**Propan-2-ol:**

Mobility : Mobile in soils

**Propan-1-ol:**

Mobility : Mobile in soils

**Biphenyl-2-ol:**

Mobility : No data available



**kodan® tincture forte colourless**

**No Change Service!**

Version 01.00  
Revision Date: 06.09.2016

Date of last issue: -  
Date of first issue: 06.09.2016

---

**12.5 Results of PBT and vPvB assessment**

**Product:**

Assessment : This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

**12.6 Other adverse effects**

**Product:**

Additional ecological information : none

**SECTION 13: Disposal considerations**

**13.1 Waste treatment methods**

Product : Dispose of the product according to the defined EWC (European Waste Code) No.  
Contaminated packaging : Take empty packaging to the recycling plant.  
Waste key for the unused product : EWC 070604  
Waste key for the unused product(Group) : Waste material of HZVA from fats, lubricants, soaps, detergents, disinfectants and personal protection products.

**SECTION 14: Transport information**

**14.1 UN number**

ADR : UN 1987  
IMDG : UN 1987  
IATA : UN 1987

**14.2 UN proper shipping name**

ADR : ALCOHOLS, N.O.S.  
(Propan-2-ol, Propan-1-ol)  
IMDG : ALCOHOLS, N.O.S.  
(Propan-2-ol, Propan-1-ol)  
IATA : Alcohols, n.o.s.  
(Propan-2-ol, Propan-1-ol)

**14.3 Transport hazard class(es)**

ADR : 3  
IMDG : 3  
IATA : 3

**14.4 Packing group**

ADR  
Packing group : III  
Classification Code : F1  
Hazard Identification Number : 30

**kodan® tincture forte colourless**

**No Change Service!**

Version                      Revision Date:  
01.00                        06.09.2016

Date of last issue: -  
Date of first issue: 06.09.2016

---

Labels                        : 3  
Tunnel restriction code    : D/E

**IMDG**

Packing group                : III  
Labels                        : 3  
EmS Code                     : F-E, S-D

**IATA**

Packing instruction (cargo aircraft) : 366  
Packing group                : III  
Labels                        : Flammable Liquid

**14.5 Environmental hazards**

**ADR**

Environmentally hazardous    : no

**IMDG**

Marine pollutant                : no

**14.6 Special precautions for user**

For personal protection see section 8.

**14.7 Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code**

Not applicable for product as supplied.

**SECTION 15: Regulatory information**

**15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**

REACH - Candidate List of Substances of Very High Concern for Authorisation (Article 59)                      : Not applicable

Regulation (EC) No 850/2004 on persistent organic pollutants                      : Not applicable

Legislation on the control of major-accident hazards involving dangerous substances                      : The product belongs to at least one of the categories 1 through 11 mentioned in Annex 1 of the Directive 1996/82/EC concerning the control of major accident hazards.

Volatile organic compounds                      : Volatile organic compounds (VOC) content: 55 %, Directive 2010/75/EC on the limitation of emissions of volatile organic compounds

Other regulations                      : Take note of Directive 98/24/EC on the protection of the health and safety of workers from the risks related to chemical agents at work. Take note of Directive 2000/39/EC establishing a first list of indicative occupational exposure limit values.

**15.2 Chemical safety assessment**

Exempt

**SECTION 16: Other information**

---

**kodan® tincture forte colourless****No Change Service!**Version                      Revision Date:  
01.00                        06.09.2016Date of last issue: -  
Date of first issue: 06.09.2016**Full text of H-Statements**

|      |   |
|------|---|
| H225 | : Highly flammable liquid and vapour.                   |
| H318 | : Causes serious eye damage.                            |
| H319 | : Causes serious eye irritation.                        |
| H335 | : May cause respiratory irritation.                     |
| H336 | : May cause drowsiness or dizziness.                    |
| H400 | : Very toxic to aquatic life.                           |
| H410 | : Very toxic to aquatic life with long lasting effects. |

**Full text of other abbreviations**

|                 |  |
|-----------------|--|
| Aquatic Acute   | : Acute aquatic toxicity                           |
| Aquatic Chronic | : Chronic aquatic toxicity                         |
| Eye Dam.        | : Serious eye damage                               |
| Eye Irrit.      | : Eye irritation                                   |
| Flam. Liq.      | : Flammable liquids                                |
| STOT SE         | : Specific target organ toxicity - single exposure |

ADN - European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways; ADR - European Agreement concerning the International Carriage of Dangerous Goods by Road; AICS - Australian Inventory of Chemical Substances; ASTM - American Society for the Testing of Materials; bw - Body weight; CLP - Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008; CMR - Carcinogen, Mutagen or Reproductive Toxicant; DIN - Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECHA - European Chemicals Agency; EC-Number - European Community number; ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. - Not Otherwise Specified; NO(A)EC - No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; RID - Regulations concerning the International Carriage of Dangerous Goods by Rail; SADT - Self-Accelerating Decomposition Temperature; SDS - Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TRGS - Technical Rule for Hazardous Substances; TSCA - Toxic Substances Control Act (United States); UN - United Nations; vPvB - Very Persistent and Very Bioaccumulative

**Further information**

Classification and procedure used to derive the classification for mixtures according to Regulation (EC) No. 1272/2008

|                    |                          |
|--------------------|--------------------------|
| Flam. Liq. 3, H226 | : On basis of test data. |
| Eye Irrit. 2, H319 | : On basis of test data. |

## Product information

on the basis of Regulation (EC) No. 1907/2006



### **kodan® tincture forte colourless**

**No Change Service!**

Version  
01.00

Revision Date:  
06.09.2016

Date of last issue: -  
Date of first issue: 06.09.2016

---

STOT SE 3, H336 : Calculation method

|| Changes compared with the previous edition!!!

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.