



Sterile ready to use solution based on non-ionic surfactants for the cleaning of surfaces in cleanrooms.

perform® sterile cleaner ND

Our Plus

- ready to use
- sterile (γ -irradiated) and microbially filtered (0.2 μm)
- available in convenient sterile 500 ml spray system
- colour and perfume free
- excellent cleaning performance

Application areas

Cleaning of surfaces and equipment / devices in cleanrooms.

Especially for Grade A and B clean areas, according to:

- EU Guidelines to Good Manufacturing Practice for medicinal products, Annex 1: Manufacture of sterile medicinal products, March 2009.
- FDA Guidance for Industry: Sterile drug products produced by aseptic processing - current Good Manufacturing Practice, U.S. Department of Health and Human Services; Food and Drug Administration, September 2004.

Instructions for use

Apply perform® sterile cleaner ND undiluted to objects and surfaces and wipe if necessary.



perform[®] sterile cleaner ND

Product data

Labelling according to Regulation (EC) No. 648/2004: < 5 % non-ionic surfactants.

Special advice

For all of the related safety information please refer to the Safety Data Sheet. To request the Safety Data Sheet please

E-mail: tech.uk@schuelke.com

The sterile 500 ml spray system offers the following advantages:

- convenient and space saving
- ergonomically shaped spray pump will not slip out of a gloved hand
- still sterile after opening

The bottles are sterile (with radiation indicator) and double-wrapped individually.

Information for order

Item	Delivery form	Item no.
perform sterile cleaner ND 500 ml FL	10/Carton	182002

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining out high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature/reports on all schulke products online at <http://www.schuelke.com/>

For any queries you might have about the product or using it please call 0114 2543500 or E-mail mail.uk@schuelke.com



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

schülke UK
Schülke & Mayr UK Ltd
Cygnet House, Meadowhall
Sheffield, S9 1AT
Phone +44-114-254 35 00
Fax +44-114-254 35 01
mail.uk@schuelke.com

schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 (0) 40 - 52100 - 0
Fax +49 (0) 40 - 52100 - 318
www.schuelke.com
mail@schuelke.com



A company of the Air Liquide-Group