According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® sterile alcohol IPA (WFI)

No Change Service!

Date of last issue: 24.01.2021

Version Revision Date: 04.00 15.09.2022

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

Trade name : perform® sterile alcohol IPA (WFI)

Unique Formula Identifier

(UFI)

W1D0-10F5-M004-FCE4

1.2 Relevant identified uses of the substance or mixture and uses advised against

Use of the Sub- : Disinfectants and general biocidal products

stance/Mixture

Recommended restrictions

on use

Diominostante ana general biocidal product

: For professional users only.

1.3 Details of the supplier of the safety data sheet

Producer : Schülke & Mayr GmbH

Robert-Koch-Str. 2

22851 Norderstedt

Germany

Telephone: +49 (0)40/ 52100-0 Telefax: +49 (0)40/ 52100318

mail@schuelke.com www.schuelke.com

Supplier : Schülke & Mayr UK Ltd.

Cygnet House

1, Jenkin Road, Meadowhall

Sheffield S9 1AT United Kingdom

Telephone: +44 114 254 35 00 Telefax: +44 114 254 35 01 mail.uk@schulke.com

E-mail address of person

responsible for the SDS/Contact person

Application Specialists +49 (0)40/ 521 00 666 AD@schuelke.com

(Schülke & Mayr UK Ltd.: +44-1142543500)

1.4 Emergency telephone number

Emergency telephone num-

ber

Carechem 24 International:+44 1235 239670

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® sterile alcohol IPA (WFI)

No Change Service!

Version Revision Date: 04.00 15.09.2022

Date of last issue: 24.01.2021

#### **SECTION 2: Hazards identification**

#### 2.1 Classification of the substance or mixture

Classification (REGULATION (EC) No 1272/2008) as amended by The Chemicals (Health and Safety) and Genetically Modified Organisms (Contained Use) (Amendment etc.) (EU Exit) Regulations 2019)

Flammable liquids, Category 2 H225: Highly flammable liquid and vapour.

Eye irritation, Category 2 H319: Causes serious eye irritation.

Specific target organ toxicity - single exposure, Category 3, Central nervous

system

H336: May cause drowsiness or dizziness.

#### 2.2 Label elements

Labelling (REGULATION (EC) No 1272/2008) as amended by The Chemicals (Health and Safety) and Genetically Modified Organisms (Contained Use) (Amendment etc.) (EU Exit) Regulations 2019)

Hazard pictograms :





Signal word : Danger

Hazard statements : H225 Highly flammable liquid and vapour.

H319 Causes serious eve irritation.

H336 May cause drowsiness or dizziness.

Supplemental Hazard

Statements

EUH066

Repeated exposure may cause skin dryness or

cracking.

Precautionary statements : Prevention:

P210 Keep away from heat, hot surfaces, sparks, open

flames and other ignition sources. No smoking.

P233 Keep container tightly closed.P261 Avoid breathing vapours/ spray.P280 Wear eye protection/ face protection.

Storage:

P403 + P235 Store in a well-ventilated place. Keep cool.

Disposal:

P501 Dispose of contents/ container to an approved waste

disposal plant.

Hazardous components which must be listed on the label:

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® sterile alcohol IPA (WFI)No Change Service!VersionRevision Date:Date of last issue: 24.01.2021

04.00 15.09.2022

propan-2-ol

#### 2.3 Other hazards

This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

Take precautionary measures against static discharges.

Vapours are heavier than air and may spread along floors.

# **SECTION 3: Composition/information on ingredients**

#### 3.2 Mixtures

Chemical nature : Alcoholic solution

### Components

Chemical name	CAS-No.	Classification	Concentration
	EC-No.		(% w/w)
	Index-No.		
	Registration number		
propan-2-ol	67-63-0	Flam. Liq. 2; H225	>= 50 - < 70
	200-661-7	Eye Irrit. 2; H319	
	603-117-00-0	STOT SE 3; H336	
	01-2119457558-25-	(Central nervous	
	XXXX	system)	

For explanation of abbreviations see section 16.

# **SECTION 4: First aid measures**

# 4.1 Description of first aid measures

General advice : Take off all contaminated clothing immediately.

Call a POISON CENTER or doctor/ physician.

If inhaled : IF INHALED: Remove person to fresh air and keep comforta-

ble for breathing.

If symptoms persist, call a physician.

In case of skin contact : Wash off with plenty of water.

If skin irritation occurs: Get medical advice/ attention.

In case of eye contact : IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue

rinsing.

If eye irritation persists: Get medical advice/ attention.

If swallowed : IF SWALLOWED: Call a POISON CENTER or doctor/ physi-

cian if you feel unwell.

## 4.2 Most important symptoms and effects, both acute and delayed

Symptoms : Treat symptomatically.

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® sterile alcohol IPA (WFI) No Change Service! **Revision Date:** Date of last issue: 24.01.2021 Version

04.00 15.09.2022

> Risks Causes serious eye irritation.

> > May cause drowsiness or dizziness.

Repeated exposure may cause skin dryness or cracking.

4.3 Indication of any immediate medical attention and special treatment needed

**Treatment** For specialist advice physicians should contact the Poisons

Information Service.

**SECTION 5: Firefighting measures** 

5.1 Extinguishing media

Suitable extinguishing media Water spray jet

Alcohol-resistant foam

Dry powder

Carbon dioxide (CO2)

Unsuitable extinguishing

media

Do NOT use water jet.

Foam

5.2 Special hazards arising from the substance or mixture

Specific hazards during fire-

fighting

Cool closed containers exposed to fire with water spray.

ucts

Hazardous combustion prod- : No hazardous combustion products are known

5.3 Advice for firefighters

for firefighters

Special protective equipment : In the event of fire, wear self-contained breathing apparatus.

**SECTION 6: Accidental release measures** 

6.1 Personal precautions, protective equipment and emergency procedures

Personal precautions Ensure adequate ventilation.

Remove all sources of ignition.

6.2 Environmental precautions

Environmental precautions Avoid subsoil penetration.

6.3 Methods and material for containment and cleaning up

Methods for cleaning up Wipe up with absorbent material (e.g. cloth, fleece).

Soak up with inert absorbent material (e.g. sand, silica gel,

acid binder, universal binder, sawdust).

6.4 Reference to other sections

see Section 8 + 13

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® sterile alcohol IPA (WFI)

No Change Service!

Date of last issue: 24.01.2021

Version Revision Date: 04.00 15.09.2022

**SECTION 7: Handling and storage** 

7.1 Precautions for safe handling

Advice on safe handling : Provide sufficient air exchange and/or exhaust in work rooms.

Advice on protection against

fire and explosion

Vapours may form explosive mixtures with air. Take measures

to prevent the build up of electrostatic charge.

Hygiene measures : Take off all contaminated clothing immediately. Keep away

from food and drink.

7.2 Conditions for safe storage, including any incompatibilities

Requirements for storage areas and containers

: Store in original container. Keep containers tightly closed in a

cool, well-ventilated place.

Further information on stor-

age conditions

Keep away from direct sunlight. Keep at temperature not ex-

ceeding 25 °C.

Advice on common storage

Do not store together with explosives, oxidizing agents, organ-

ic peroxides and infectious products.

7.3 Specific end use(s)

Specific use(s) : none

# **SECTION 8: Exposure controls/personal protection**

### 8.1 Control parameters

#### **Occupational Exposure Limits**

Components	CAS-No.	Value type (Form of exposure)	Control parameters	Basis
propan-2-ol	67-63-0	TWA	400 ppm 999 mg/m3	GB EH40
		STEL	500 ppm 1,250 mg/m3	GB EH40

#### **Derived No Effect Level (DNEL):**

Substance name	End Use	Exposure routes	Potential health effects	Value
propan-2-ol	Workers	Skin contact	Long-term systemic effects	888 mg/kg
	Workers	Inhalation	Long-term systemic effects	500 mg/m3

### **Predicted No Effect Concentration (PNEC):**

Substance name	Environmental Compartment	Value
propan-2-ol	Fresh water	140.9 mg/l
	Marine water	140.9 mg/l
	Fresh water sediment	552 mg/kg
	Marine sediment	552 mg/kg

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® sterile alcohol IPA (WFI) No Change Service!

Version Revision Date: Date of last issue: 24.01.2021

04.00 15.09.2022

II	Soil	28 mg/kg
	Intermittent use/release	140.9 mg/l
	Effects on waste water treatment plants	2251 mg/l
	Oral	160 mg/kg food

#### 8.2 Exposure controls

Personal protective equipment

Eye/face protection : Safety glasses with side-shields conforming to EN166

Hand protection

Directive : The selected protective gloves have to satisfy the specifica-

tions of Regulation (EU) 2016/425 and the standard EN 374

derived from it.

Remarks : Splash protection: disposable nitrile rubber gloves e.g.

Dermatril (layer thickness: 0.11 mm) made by KCL or gloves from other manufacturers offering the same protection. Prolonged contact: Nitrile rubber gloves e.g. Camatril (>120 Min., layer thickness: 0.40 mm) or butyl rubber gloves e.g. Butoject (>480 Min., layer thickness: 0.70 mm) made by KCL or gloves from other manufacturers offering the same protec-

tion.

Skin and body protection : Work uniform or laboratory coat.

Respiratory protection : If the occupational exposure limits cannot be met, in excep-

tional cases suitable respiratory equipment should be worn

only for a short period of time. Recommended Filter type:

A-P2 or ABEK-P2

Respiratory protection complying with EN 141.

Protective measures : Avoid contact with skin and eyes.

Do not breathe vapours or spray mist.

## **SECTION 9: Physical and chemical properties**

## 9.1 Information on basic physical and chemical properties

Appearance : liquid

Colour : colourless

Odour : alcohol-like

Odour Threshold : not determined

pH : Not applicable

Melting point/freezing point : < -5 °C

Decomposition temperature Not applicable

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® sterile alcohol IPA (WFI) No Change Service!

Version Revision Date: Date of last issue: 24.01.2021

04.00 15.09.2022

Boiling point/boiling range : 82 - 83 °CMethod: DIN 53171

Raw material

Flash point : 22 °C

Method: DIN 51755 Part 1

Evaporation rate : No data available

Upper explosion limit / Upper

flammability limit

12 %(V) Raw material literature value

Lower explosion limit / Lower

flammability limit

: 2 %(V)

Raw material literature value

Vapour pressure : > 42 hPa (20 °C)

Raw material literature value

Relative vapour density : No data available

Density : 0.87 - 0.881 g/cm3 (20 °C)

Solubility(ies)

Water solubility : completely soluble (20 °C)

Partition coefficient: n-

octanol/water

Not applicable

Viscosity

Viscosity, dynamic : > 2.5 mPa\*s (20 °C)

Raw material

Explosive properties : Vapours may form explosive mixtures with air.

Oxidizing properties : The substance or mixture is not classified as oxidizing.

9.2 Other information

Flammability (liquids) : Sustains combustion

Metal corrosion rate : None reasonably foreseeable.

Self-ignition : No data available

# **SECTION 10: Stability and reactivity**

## 10.1 Reactivity

No dangerous reaction known under conditions of normal use.

# 10.2 Chemical stability

The product is chemically stable.

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® sterile alcohol IPA (WFI) No Change Service!

Version Revision Date: Date of last issue: 24.01.2021

04.00 15.09.2022

10.3 Possibility of hazardous reactions

Hazardous reactions : Vapours may form explosive mixture with air.

10.4 Conditions to avoid

Conditions to avoid : Heat, flames and sparks.

Exposure to sunlight.

10.5 Incompatible materials

Materials to avoid : Strong acids and oxidizing agents

10.6 Hazardous decomposition products

No decomposition if stored and applied as directed.

**SECTION 11: Toxicological information** 

11.1 Information on toxicological effects

**Acute toxicity** 

Not classified based on available information.

**Components:** 

propan-2-ol:

Acute oral toxicity : LD50 (Rat): 5,840 mg/kg

Acute inhalation toxicity : LC50 (Rat): 39 mg/l

Exposure time: 4 h
Test atmosphere: vapour

Acute dermal toxicity : LD50 (Rabbit): 13,900 mg/kg

Method: OECD Test Guideline 402

Skin corrosion/irritation

Repeated exposure may cause skin dryness or cracking.

**Components:** 

propan-2-ol:

Result : No skin irritation

Serious eye damage/eye irritation

Causes serious eye irritation.

Components:

propan-2-ol:

Result : Eye irritation

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® sterile alcohol IPA (WFI)No Change Service!VersionRevision Date:Date of last issue: 24.01.2021

04.00 15.09.2022

## Respiratory or skin sensitisation

#### Skin sensitisation

Not classified based on available information.

#### Respiratory sensitisation

Not classified based on available information.

#### **Components:**

### propan-2-ol:

Test Type : Buehler Test Species : Guinea pig

Result : Did not cause sensitisation on laboratory animals.

### Germ cell mutagenicity

Not classified based on available information.

### **Components:**

### propan-2-ol:

Genotoxicity in vitro : Test Type: Ames test

Method: Mutagenicity (Escherichia coli - reverse mutation

assay)

Result: Non mutagenic

Genotoxicity in vivo : Species: Mouse

Method: Mutagenicity (micronucleus test)

Result: Non mutagenic

Germ cell mutagenicity- As-

sessment

: Not mutagenic in Ames Test

#### Carcinogenicity

Not classified based on available information.

### **Components:**

#### propan-2-ol:

Remarks : Based on available data, the classification criteria are not met.

### Reproductive toxicity

Not classified based on available information.

#### **Components:**

# propan-2-ol:

Effects on foetal develop- : Species: Rat

ment Application Route: Oral

General Toxicity Maternal: NOAEL: 400 mg/kg body weight

Reproductive toxicity - As-

sessment

: Based on available data, the classification criteria are not met.

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® sterile alcohol IPA (WFI) No Change Service!

**Revision Date:** Date of last issue: 24.01.2021

04.00 15.09.2022

STOT - single exposure

May cause drowsiness or dizziness.

Components:

propan-2-ol:

Assessment May cause drowsiness or dizziness.

STOT - repeated exposure

Not classified based on available information.

Components:

propan-2-ol:

Remarks Based on available data, the classification criteria are not met.

Repeated dose toxicity

**Components:** 

propan-2-ol:

Remarks No data available

**Aspiration toxicity** 

Not classified based on available information.

**Further information** 

**Product:** 

Remarks No data is available on the product itself.

### **SECTION 12: Ecological information**

## 12.1 Toxicity

#### Components:

propan-2-ol:

Toxicity to fish LC50 (Pimephales promelas (fathead minnow)): 9,640 mg/l

Exposure time: 96 h

aquatic invertebrates

Toxicity to daphnia and other : EC50 (Daphnia magna (Water flea)): 10,000 mg/l

Exposure time: 48 h

Toxicity to algae/aquatic

plants

: EC50 (Desmodesmus subspicatus (green algae)): > 100 mg/l

Exposure time: 72 h

Test Type: static test

EC50 (green algae): 1,800 mg/l

Exposure time: 7 d

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® sterile alcohol IPA (WFI) No Change Service!

**Revision Date:** 

Date of last issue: 24.01.2021

04.00 15.09.2022

## 12.2 Persistence and degradability

**Product:** 

Version

Biodegradability : Result: Readily biodegradable.

**Components:** 

propan-2-ol:

Biodegradability : Result: Readily biodegradable.

12.3 Bioaccumulative potential

**Components:** 

propan-2-ol:

Bioaccumulation : Remarks: No bioaccumulation is to be expected (log Pow <=

4).

Partition coefficient: n-

: log Pow: 0.05 (20 °C)

octanol/water

Method: OECD Test Guideline 107

12.4 Mobility in soil

**Components:** 

propan-2-ol:

Mobility : Remarks: Mobile in soils

12.5 Results of PBT and vPvB assessment

**Product:** 

Assessment : This substance/mixture contains no components considered

to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of

0.1% or higher.

12.6 Other adverse effects

**Product:** 

Endocrine disrupting poten-

tial

The substance/mixture does not contain components considered to have endocrine disrupting properties according to

REACH Article 57(f) or Commission Delegated regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605 at

levels of 0.1% or higher.

Additional ecological infor-

mation

: No data is available on the product itself.

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® sterile alcohol IPA (WFI) No Change Service!

Version Revision Date: Date of last issue: 24.01.2021

04.00 15.09.2022

## **SECTION 13: Disposal considerations**

#### 13.1 Waste treatment methods

Product : Disposal together with normal waste is not allowed. Special

disposal required according to local regulations.

Contaminated packaging : Empty containers should be taken to an approved waste han-

dling site for recycling or disposal.

# **SECTION 14: Transport information**

# 14.1 UN number

 ADR
 : UN 1219

 IMDG
 : UN 1219

 IATA
 : UN 1219

14.2 UN proper shipping name

ADR : ISOPROPANOL

IMDG : ISOPROPANOL, SOLUTION

IATA : Isopropanol, solution

### 14.3 Transport hazard class(es)

 ADR
 : 3

 IMDG
 : 3

 IATA
 : 3

### 14.4 Packing group

## ADR

Packing group : II
Classification Code : F1
Hazard Identification Number : 33
Labels : 3
Tunnel restriction code : (D/E)

**IMDG** 

Packing group : II
Labels : 3
EmS Code : F-E, S-D

IATA (Cargo)

Packing instruction (cargo : 364

aircraft)

Packing instruction (LQ) : Y341
Packing group : II

Labels : Flammable liquid

IATA (Passenger)

Packing instruction (passen: 353

ger aircraft)

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



 perform® sterile alcohol IPA (WFI)
 No Change Service!

 Version
 Revision Date:

 Date of last issue: 24.01.2021

04.00 15.09.2022

Packing instruction (LQ) : Y341
Packing group : II

Labels : Flammable liquid

14.5 Environmental hazards

ADR

Environmentally hazardous : no

**IMDG** 

Marine pollutant : no

14.6 Special precautions for user

The transport classification(s) provided herein are for informational purposes only, and solely based upon the properties of the unpackaged material as it is described within this Safety Data Sheet. Transportation classifications may vary by mode of transportation, package sizes, and variations in regional or country regulations.

14.7 Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable for product as supplied.

**SECTION 15: Regulatory information** 

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

Relevant EU provisions transposed through retained EU law

UK REACH List of restrictions (Annex 17) : Conditions of restriction for the fol-

lowing entries should be considered:

Number on list 3 Not applicable

UK REACH Candidate list of substances of very high

concern (SVHC) for Authorisation

The Persistent Organic Pollutants Regulations (retained

Regulation (EU) 2019/1021 as amended for Great Brit-

ain)

: Not applicable

Regulation (EC) No 1005/2009 on substances that de-

plete the ozone layer

Not applicable

UK REACH List of substances subject to authorisation

(Annex XIV)

Not applicable

Registration number : EU-0023656-0000

Volatile organic compounds : Directive 2010/75/EU of 24 November 2010 on industrial

emissions (integrated pollution prevention and control) Volatile organic compounds (VOC) content: 63.1 %

The components of this product are reported in the following inventories:

TCSI : On the inventory, or in compliance with the inventory

TSCA : All substances listed as active on the TSCA inventory

AIIC : On the inventory, or in compliance with the inventory

Z11211 ZSDB\_P\_GB EN

Page 13/15

According to REACH etc. (Amendment etc.) (EU Exit) Regulations



perform® sterile alcohol IPA (WFI) **Revision Date:** 

No Change Service! Date of last issue: 24.01.2021

04.00 15.09.2022

Version

DSL All components of this product are on the Canadian DSL

**ENCS** On the inventory, or in compliance with the inventory

**ISHL** On the inventory, or in compliance with the inventory

**KECI** On the inventory, or in compliance with the inventory

**PICCS** On the inventory, or in compliance with the inventory

**IECSC** On the inventory, or in compliance with the inventory

**NZIoC** Not in compliance with the inventory

**TECI** On the inventory, or in compliance with the inventory

#### 15.2 Chemical safety assessment

Exempt

#### **SECTION 16: Other information**

## **Full text of H-Statements**

Highly flammable liquid and vapour. H319 Causes serious eye irritation. H336 May cause drowsiness or dizziness.

#### Full text of other abbreviations

Eve Irrit. Eve irritation Flam. Liq. Flammable liquids

Specific target organ toxicity - single exposure STOT SE UK. EH40 WEL - Workplace Exposure Limits GB EH40

GB EH40 / TWA Long-term exposure limit (8-hour TWA reference period) GB EH40 / STEL Short-term exposure limit (15-minute reference period)

ADN - European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways; ADR - Agreement concerning the International Carriage of Dangerous Goods by Road; AIIC - Australian Inventory of Industrial Chemicals; ASTM - American Society for the Testing of Materials; bw - Body weight; CLP - Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008; CMR - Carcinogen, Mutagen or Reproductive Toxicant; DIN - Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECHA -European Chemicals Agency; EC-Number - European Community number; ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - Interna-

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® sterile alcohol IPA (WFI)

No Change Service!

Version Revision Date: 04.00 15.09.2022

Date of last issue: 24.01.2021

tional Convention for the Prevention of Pollution from Ships; n.o.s. - Not Otherwise Specified; NO(A)EC - No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; RID - Regulations concerning the International Carriage of Dangerous Goods by Rail; SADT - Self-Accelerating Decomposition Temperature; SDS - Safety Data Sheet; SVHC - Substance of very high concern; TCSI - Taiwan Chemical Substance Inventory; TECI - Thailand Existing Chemicals Inventory; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative

#### **Further information**

#### Classification of the mixture:

Classification procedure:

Flam. Liq. 2 H225 Based on product data or assessment

Eye Irrit. 2 H319 Calculation method STOT SE 3 H336 Calculation method

Changes since the last version are highlighted in the margin. This version replaces all previous versions.

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.