



Spray bottle holder for clean rooms.

perform[®] spray bottle holder

Our plus:

- Compatible with all spray bottles from the perform[®] range
- Fully autoclavable
- Made of high quality stainless steel
- Easy installation

Function description

The pharma spray bottle holder can be specifically used for glove disinfection in conjunction with the products of the perform[®] range. The dispenser is compatible with the 500 ml and 1 L ready-to-use spray bottles for example from the products perform[®] sterile alcohol or IPA perform[®] sterile alcohol EP.

The pharma spray bottle holder is made of stainless steel and is autoclavable.

Area of application

The pharma spray bottle holder provides an easy to use application method that is particularly suitable for Grade A / B clean-rooms.

Installation

There are two ways to install the pharma spray bottle holder:

- The bottle holder can be mounted directly via the two screw holes.
- The pharma spray bottle holder mounting plate is glued to the wall and the bottle holder itself is attached via the two pins on the mounting plate. The bottle holder can thus be easily removed for cleaning and disinfection. For attaching the mounting plate we recommend the glabete adhesive system.

Mounting via adhesive and plate facilitates the regular autoclaving of the pharma spray bottle holder

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining out high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: www.schuelke.com. For individual questions:

Customer Sales Service
Phone: +49 40 52100-666
E-Mail: info@schuelke.com



Schülke & Mayr GmbH is certified according to DIN EN ISO 9001, DIN EN ISO 14001 and DIN EN ISO 13485 (Reg.-No. 004567-MP23) and has a validated environmental system in accordance with the Eco Audit Regulation (Reg.-No. DE-150-00003).

Schülke & Mayr GmbH
22840 Norderstedt, Germany
Phone +49 (0) 40 - 521 00 - 0
Fax +49 (0) 40 - 521 00 - 318
www.schuelke.com
mail@schuelke.com