



Flexible bottle holder for easy mounting on tube frames.

Point of Care Halter

Our Plus

- highly flexible
- can be mounted on horizontal and vertical tubular frames
- nearly stepless assembly
- with anti-slip surface for optimal grip
- stable and tear-resistant material
- can be thermally reprocessed up to 60°

Application area

- ensures the availability of hand rubs directly at the Point of Care
- suitable for 500ml-hyclick® bottles and 500ml-bottles with separate screw-on pump
- can be quickly mounted on horizontal and vertical tubular frames (e.g. patient bed) with a diameter of Ø 18 - 40 mm
- can be wiped down with a soft disinfection wipe (e.g. with mikrozid® universal wipes green line) or thermally reprocessed up to 60°C (e.g. with thermosept® X-tra enzymatic cleaner)

Information for order

Item	Delivery form	Item no.
PoC Holder Manual DE/EN	1	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

Related products

- **desderman®**
- **desmanol® care**

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

schülke Headquarters
Schülke & Mayr GmbH
 Robert-Koch-Straße 2
 22851 Norderstedt | Germany
 Phone +49 40 52100582
www.schuelke.com
info@schuelke.com



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.