



Neutraliser based on citric acid for the automated reprocessing of medical devices and laboratory accessories.

thermosept® NKZ

Our Plus

- phosphate free
- good material compatibility
- optimises cleaning results after neutral cleaning
- prevents the build-up of deposits

Application areas

thermosept® NKZ is used after alkaline cleaning for the treatment of surgical, microsurgery and minimal invasive surgery instruments, stainless steel containers, rigid endoscopes, anaesthesia equipment, operating theatre shoes and baby bottles, in order to prevent the possibility of alkaline residues. Furthermore it is suitable for use after neutral cleaning with thermosept® RKN-zym for improved cleaning results.

thermosept® NKZ can be used in cycle assembly lines and individual machines.

Instructions for use

After neutral cleaning thermosept® NKZ can be used at a concentration of 0.1% to optimise the results for cleaning. The use of demineralised water during the final rinse will improve rinsing results.



Product data

Composition:
thermosept NKZ contains organic acids. The product contains no surfactants or phosphates.

Special advice

For all of the related safety information please refer to the Safety Data Sheet. To request the Safety Data Sheet please E-mail: tech.uk@schuelke.com
Do not swallow.

Information for order

Item	Delivery form	Item no.
thermosept NKZ 5 kg KA	1/Canister	146503

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature/reports on all schulke products online at <http://www.schuelke.com/>

For any queries you might have about the product or using it please call 0114 2543500 or E-mail mail.uk@schuelke.com



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.



schülke UK
Schülke & Mayr UK Ltd
Cygnet House, Meadowhall
Sheffield, S9 1AT
Phone +44-114-254 35 00
Fax +44-114-254 35 01
mail.uk@schuelke.com

schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com