



Unique alcohol formulation for virucidal protection for the cleaning and disinfection of hard surfaces

## mikrozid<sup>®</sup> AF wipes premium

### Our Plus

- blend of two alcohols for a very high performance against viruses within 1 min contact time
- cleaning and disinfection, 2 complete actions with 1 product
- added surfactants for boosted cleaning performance
- broad spectrum of efficacy
- high quality fleece wipes for double the coverage\*

### Application areas

Alcohol-based rapid disinfection of medical devices and surfaces in all areas where there is an increased risk of infection and where short contact times are required, e.g.:

- Patient treatment units
- Operating tables and nearby work surfaces
- Surfaces of medical equipment

### Instructions for use

Wipe areas thoroughly with the impregnated cloth; ensure thorough wetting and leave to work. If necessary, use several wipes for large areas. Make sure that all visible dirt is removed prior to disinfection. Make sure to close the lid tightly after use. Do not treat surfaces that are sensitive to alcohol, e.g. acrylic glass. Not suitable for the disinfection of invasive medical devices. schülke guarantees a shelf life of 3 months after the package has been opened.

### Microbiological efficacy

Efficacy	Concentration	Contact time
<b>bactericidal</b> <i>in accordance with VAH Guidelines with short contact times</i>	ready-to-use	1 min.
<b>MRSA</b>	ready-to-use	1 min.

Efficacy	Concentration	Contact time
<b>tuberculocidal</b> <i>EN14348 - dirty conditions</i>	ready-to-use	1 min.
<b>yeastcidal</b> <i>in accordance with VAH Guidelines with short contact times</i>	ready-to-use	1 min.
<b>virucidal against enveloped viruses</b> <i>- dirty conditions</i>	ready-to-use	1 min.
<b>Vaccinia virus</b>	ready-to-use	30 sec.
<b>Norovirus</b> <i>EN14476 - dirty conditions</i>	ready-to-use	1 min.
<b>Rotavirus</b>	ready-to-use	30 sec.

CE 0297



# mikrozid<sup>®</sup> AF wipes

premium

## Product data

Composition:

100 g solution contains the following active substances:  
25 g Ethanol (94 %), 35 g Propan-1-ol

## Special advice

For all of the related safety information please refer to the Safety Data Sheet. To view the Safety Data Sheet please visit [www.schuelke.com](http://www.schuelke.com)

**Use disinfectant safely. Always read the label and product information before use.**

The product is highly compatible with metals and synthetics (except acrylic glass and alcohol-sensitive varnishes). Please follow carefully any national fire and explosion protection guidelines when using alcohol-based disinfectants. Only treat small surfaces. Surface coverage\* approximately 0.5m<sup>2</sup>. The wiping range depends on the room temperature and on the structure of the surface to be disinfected. More information available on request.

## Information for order

Item	Delivery form	Item no.
mikrozid AF wipes prem 50 ST SP	12/Carton	70000808

## Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

## Expert opinion and information

Please visit our website for an overview of all available literature/reports on all schulke products online at <http://www.schuelke.com/>

For any queries you might have about the product or using it please call 01 14 2543500 or E-mail [mail.uk@schuelke.com](mailto:mail.uk@schuelke.com)

Item no.	GTIN no.
70000808	4032651959075

\*Comparison is based on standard mikrozid<sup>®</sup> wipes. The surface coverage of one wipe depends on the surrounding conditions and the structure of the surface.



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

**schülke UK**  
Schülke & Mayr UK Ltd  
Cygnets House, Meadowhall  
Sheffield, S9 1AT  
Phone +44-114-254 35 00  
Fax +44-114-254 35 01  
[mail.uk@schuelke.com](mailto:mail.uk@schuelke.com)

**schülke Headquarters**  
Schülke & Mayr GmbH  
Robert-Koch-Str. 2  
22851 Norderstedt, Germany  
Phone +49 40 - 52100 - 0  
Fax +49 40 - 52100 - 318  
[www.schuelke.com](http://www.schuelke.com)  
[info@schuelke.com](mailto:info@schuelke.com)