



Examination glove made from an advanced material that is inner coated with an aloe and vitamin E conditioning lotion

Sempercare® Nitrile Aloe Gloves

Our Plus

- contains aloe to improve skin moisture
- contains vitamin E for improved skin health
- excellent tactile sensitivity due to the reduced wall thickness
- safe alternative for latex type I allergy sufferers
- soft and elastic advanced nitrile material
- lavender blue colour

Product Profile

A hand-conditioning lotion containing aloe and vitamin E is applied to the inside of the gloves helping to improve skin moisture and health. Aloe can nourish your skin in a natural way, which makes it ideal for dry and chapped skin. The structured finger tips guarantee a sure grip, even in a moist environment. The lavender-blue colour helps distinguish the gloves from latex alternatives.

The gloves meet European impermeability standards with proven protection from the transmission of infectious microorganisms. As a latex free alternative, the gloves provide a safe alternative for latex type I allergy sufferers and are low in residual chemicals.

CE 0121



Sempercare[®] Nitrile Aloe Gloves

Special advice

For all of the related safety information please refer to the Safety Data Sheet. To request the Safety Data Sheet please E-mail: tech.uk@schuelke.com

Examination glove and disposable protective gloves made of a nitrile butadiene rubber (NBR), powder-free, non sterile. Ambidextrous. Free from natural latex Synthetic copolymer. Store in original packaging in a dry and dark place at 10°C to 30°C. Protect from ozone. Dispose of in accordance with valid regulations.

This product contains Aloe Vera - not to be used in the case of hypersensitivity to this substances. For single use. Check for damage before use, do not use damaged gloves.

Information for order

Item	Delivery form	Item no.
Sempercare Nitrile Aloe - Sm 200 ST BX	10/Carton	352066
Sempercare Nitrile Aloe - XS 200 ST BX	10/Carton	352065
Sempercare Nitrile Aloe - L 200 ST BX	10/Carton	352068
Sempercare Nitrile Aloe - Med 200 ST BX	10/Carton	352067
Sempercare Nitrile Aloe - XL 180 ST BX	10/Carton	352069

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature / reports on all schülke products online at <http://www.schuelke.com/>

For any queries you might have about the product or using it please call on 0114 2543500 or E-mail mail.uk@schuelke.com



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.



schülke UK
Schülke & Mayr UK Ltd
Cygnat House, Meadowhall
Sheffield, S9 1AT
Phone +44-114-254 35 00
Fax +44-114-254 35 01
mail.uk@schuelke.com

schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com